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Gear pump drive Microprocessor controlled

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Caution

Rotating pumphead creates a pinch and crush hazard.

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Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- ▶ The circuit between mains supply and pump has to be earthed.
- ▶ The pump must not be operated outside the destined operating and environmental conditions.
- ▶ Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- ▶ The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes
- ▶ The pump must be switched off when changing the pump-head.
- ▶ The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.

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Please note

Never pump media containing particles.

If necessary, filter the media before starting your application.

Never run the pump-head under dry conditions.

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Safety precautions

- ▶ Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- ▶ We advise you never pump liquid containing particles. Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- ▶ Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.
- ▶ Only new fuses, according to the specifications stated on Page 9 in this manual, must be used.
- ▶ The fuse-holder must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Pump drive: 2 years

Pump-head: 1 year

The pump-heads are covered by the specific terms of Micropump Ltd. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use
- ▶ unauthorized modification or misuse by the user or by a third-party
- ▶ improper site preparation and maintenance
- ▶ operation outside of the environmental and electrical specifications for the product
- ▶ use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- ▶ MCP-Z Standard drive with 1 magnet. Order No. ISM 405
- ▶ 1 power cord with an IEC 320 connector (female plug) and a country specific mains plug
- ▶ 2 tubing clamps
- ▶ 1 non-return valve (PP)
- ▶ Operating manual

If ordered:

- ▶ Pump-head with 2 tubing adapters

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

➡ Complaints can only be accepted within 8 days from receipt of the goods.

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

Rear panel

1. Socket for footswitch and hand dispenser
2. Ventilator
3. Fuse-holder with voltage selector 115/230 V
4. Mains socket
5. RS232 In (female)
6. RS232 Out (male)
7. Valve connector
8. Analog interface

- input for:
- speed control
0 – 5 V or 0 – 10 V, and
0 – 20 mA or 4 – 20 mA
 - Run/Stop
- Output for:
– speed 0 – 10 V_{DC}
or 0 – 12 kHz

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Mains voltage

Mains voltage	Voltage setting	Fuse rating
220–240 V _{AC} 110–120 V _{AC}	230 V 50/60Hz 115 V 50/60Hz	2 x T1.25A/250V 2 x T2.50A/250V

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

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


Voltage setting 115/230 V Changing the fuses

- ➔ Switch the pump off, pull out the mains plug.
- 1. Pull out the fuse-holder by opening it with a small screw-driver (size 0).
- 2. Take out the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.
- 3. Insert 2 new fuses
 - ➔ 230 V: 2 x T1.25A/250V
 - ➔ 115 V: 2 x T2.5A/250V

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.
- 4. Shut the fuse-holder. The voltage rating is visible in the window.

Changing the fuses on the control board

- The 2 fuses are fixed to the control board as illustrated opposite.
- a. T5A/250V
 - b. F1.6A/250V
- ➔ Make sure that the pump is disconnected from the mains supply.
 - ⚠ **The instrument should only be opened by a qualified technician!** Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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
Starting the pump

- ▶ Check the voltage rating in the window of the fuse-holder (on rear panel). If necessary, change the rating as indicated on Page 9.
 - ▶ Mount the pump-head as stated on Pages 46 – 47.
 - ▶ Before screwing the tubing adapters into the pump-head, we recommend you wrap **PTFE tape around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.
 - ▶ Connect the tubing to your system. Fasten the tubing to the adapters by means of the tubing clamps.
 - ▶ Connect the pump to the mains and switch it on with the power supply switch.
- ➡ Enter the ID Codes of the mounted pump-heads in the program memory currently used (see Pages 15/49/50).
- ⚠ **Before you start pumping, fill the pump-head with the medium or water. Never run the pump-head under dry conditions.**
- ⚠ **Never pump liquids that contain particles. If necessary, filter the medium before use.**

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Operating panel

1. Mains switch (on/off)
2. Digital LED display
3. Reduce value
4. Increase value
5. Run / Stop
6. Calibration / Data saving
7. MAX/RESET (Multi-function, see Page 12)
8. Selecting the operating mode
➔ MODE key

9. LEDs for active operating mode

- ▶ **PUMP rpm**
Pumping by speed 60 – 6000 rpm
- ▶ **PUMP Flow rate**
Pumping by flow rate
- ▶ **DISP Time**
Dispensing by time
- ▶ **DISP Volume**
Dispensing by volume
- ▶ **PAUSE Time**
Pause for intermittent dispensing
- ▶ **PROGRAM**
Program selection 1 – 4
- ▶ **TOTAL**
Read-out of totally delivered volume

➔ The drive runs only in clock-wise direction.

⚠ Never run the pump-head under dry conditions.

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

Control keys

- a. **RUN/STOP**
Starts and stops the pump.
- b. **MODE**
Changes between operating modes (see Page 11).
- c. **CAL**
Key for saving a set value.
- d. **MAX/RESET**
Multi-function key for:
- ➔ max. speed (when pump is running)
 - ➔ Reset button resets the calibration to the standard values (default, see Page 23)
 - ➔ Reset button for accumulated volume in Mode TOTAL
- e. **Increment/Decrement keys**
▲ = increase value
▼ = decrease value
Maintaining pressure on ▲ or ▼ key changes display read-out into fast mode.

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Start-up information
<p>After switching on the power supply, the following values are displayed:</p> <ol style="list-style-type: none"> 1. LED test »8.8.8.8.« 2. Firmware version 3. Baud rate 4. Identification Code (ID-Code) of the pump-head entered in the program last used (e.g.120) 5. Program last used. 6. Displays the operating mode last used. <p>⚠ Before starting the pump the first time or after mounting a new pump-head, the ID-Code must be entered <u>in each program</u> used with this pump-head (see Pages 15/49/50).</p>

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Resetting the parameters of the currently used program to the default values.

1. Power switch »OFF«
2. Push the **MODE key** while switching the drive on

The following parameters are reset:

Mode:	PUMP rpm
Speed:	1000 rpm
Dispensing time:	4.5 sec.
Dispensing volume and flow rate:	depending on pump-head
Pause time:	2.0 sec.
Number of dispensing cycles:	12

➔ **The pump-head entered in the currently set program remains stored.**

Resetting the parameters of all 4 programs to the default values

1. Power switch »OFF«
2. Push the **▼ key** while switching the drive on.

The following parameters are reset (in all 4 programs):

Pump-head: Model 120 as well as all parameters mentioned above

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Entering the pump-head identification number

In order to obtain correct pumping and dispensing values the individual identification code (ID-Code) of the mounted pump-head must be entered (**▲ for each of the 4 programs**), **especially** when a new pump is used for the first time and each time the pump-head is changed.

1. Power switch »OFF«
2. Keep CAL-key pressed and switch power »ON«
3. Press ▲ or ▼ key till the correct ID-Code of the mounted pump-head is set. The ID code is indicated on the pump-head or in the table on Page 49/50.
4. Enter the correct ID-Code with the CAL-key (MODE display switches automatically to PUMP rpm)
5. Change to PROGRAM with the MODE key and select the next program with the CAL and ▲▼ keys (confirm with CAL).


➔ Repeat procedure 1 – 4 till the ID code of the mounted pump-head is entered in all 4 programs.

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Program selection

When switching the pump on, it always selects the previously used program.

1. Change mode to PROGRAM by using the MODE-key.
2. Press CAL button, currently set program starts to blink. Change to the required program with the **▲▼** buttons.
3. Confirm with the CAL-key
4. The pump automatically returns to the last entered operating parameters of the selected program.

⚠ From now on, any changes carried out in the operating modes are automatically stored in the currently selected program.

TOTAL Volume

Change with MODE key to TOTAL. The totally delivered volume is displayed.

Three display readings are available:

Microlitre: e.g. 17.0 µl = **17.0 µ**

Millilitre: e.g. 1700 ml = **17.00**

Litre: e.g. 17.0 l = **17.0L**

(Display in litres from 9999 ml)

For zero-setting, either

- ▶ press the reset key twice, or
- ▶ switch the pump OFF.

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

Pumping by drive speed

- Change mode to PUMP rpm (60 – 6000 rpm) digitally adjustable in stops of 1 rpm
- Enter required speed with ▲▼ keys
 - ▶ Start the pump with RUN/STOP
 - ➡ The speed can also be adjusted while the pump is running.
- Maintain pressure on MAX-key while the pump is running.
 - ➡ Pump is running at maximum speed (ideal for fast filling or emptying the system).

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Pumping by flow rate

1. Change mode to PUMP Flow rate

1. Change mode to PUMP Flow rate
2. Enter the required flow rate with ▲▼ keys (is displayed in ml/min)

Start pump with RUN/STOP

- ➔ The flow rate can also be adjusted while the pump is running.
- ➔ For an accurate flow rate we recommend to calibrate the pump (see Page 19).
- ➔ During the pumping process the following values can be retrieved with the MODE-key:
 - Rotation speed (PUMP rpm)
 - Totally pumped volume (TOTAL)

3. Maintain pressure on MAX-key while the pump is running

- ➔ Pump is running at maximum speed (ideal for fast filling or emptying of the system).

⚠ Entering the ID-Code of the currently mounted pump-head (Pages 15/49/50) allows the user to set the flow rate in ml/min. This, however, is an approximate value as not yet calibrated.

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Calibrating the flow rate

- Change mode to PUMP Flow rate. Enter the required flow rate with the ▲▼ keys
- Change mode to DISP Time
 - Enter 60 seconds by using the ▲▼ keys
 - Start pump with RUN/STOP

➡ Pump stops automatically after 60 seconds

➡ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:
- Return to mode PUMP Flow rate
- Press the CAL-button (displayed value blinks) Enter the weighed or measured value with the ▲▼ keys.
- Confirm with the CAL-key (the flow rate setting returns automatically to the initially entered set point).
- Start with the RUN/STOP-key

➡ Depending on the application conditions repeated calibration may be necessary.

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
- ▶ Drehzahl (PUMP rpm)
- ▶ Fließrate (PUMP Flow rate)
- ▶ total gefördertes Volumen (TOTAL)

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Dispensing by time

The dispensing time can be entered from 0.1 s to 999 h.

1. Change mode to DISP Time

2. Enter the required dispensing time with the ▲▼ keys

➡ The display accelerates when pressure on the ▲ or ▼ key is maintained.

Seconds: 0.1" – 899.9" (in 0.1 s steps)
 Minutes: 15' – 899' (in 1 min steps)
 Hours: 15 h – 999 h (in 1 h steps)

3. Start pump with RUN/STOP

➡ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point only takes effect from the subsequent dispensing cycle.

➡ During the dispensing process the following values can be retrieved by pressing the MODE-key:

- ▶ speed (PUMP rpm)
- ▶ flow rate (PUMP Flow rate)
- ▶ totally dispensed volume (TOTAL).

➡ With the ▲▼ keys the speed and flow rate can be changed even during the dispensing process.

➡ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible.

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Dispensing by volume

1. Change mode to DISP Volume

2. Use the ▲▼ keys for entering the required dispensing volume. Three display readings are available:

Micro litre: e.g. 17.0 µl = **17.0 µ**
 Millilitre: e.g. 17.00 ml = **17.00**
 Litre: e.g. 17.0 l = **17.0L**
 (Display in litres from 9999 ml)

➔ For dispensing accurately we recommend you calibrate the pump (see Page 22)

3. Start pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

➔ The dispensing speed can be adjusted in the modes PUMP rpm and PUMP flow rate.

➔ With the ▲▼ keys the dispensing volume can also be changed even during the dispensing process. The newly entered set point appears shortly in the display taking effect, however, only with the subsequent dispensing step.

➔ During the dispensing process the following values can be retrieved with the MODE-key:

- ▶ Rotation speed (PUMP rpm)
- ▶ Flow rate (PUMP Flow rate)
- ▶ totally dispensed volume (TOTAL)

➔ With the ▲▼ keys the rotation speed and flow rate can be changed even during the dispensing process.

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PUMP	o rpm	o Flow Rate
DISP	o Time	* Volume
PAUSE	o Time	o PROGRAM
TOTAL	o	

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Calibrating the volume

- Change MODE to DISP Volume
 ➔ Use the ▲▼ keys for entering the required dispensing volume.
- Start with RUN/STOP
 ➔ Pump stops automatically.
 ➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes.
 ➔ Measure the dispensed liquid by volume or weight.
- Stay in the same mode (DISP Volume) and press the CAL-key (display blinks).
- Enter the effectively dispensed volume with the ▲▼ keys.
- Save setting with CAL-key
 ➔ Pump is calibrated and returns automatically to the initially entered set point.
- Start pump with RUN/STOP.
 ➔ Depending on the application conditions repeated calibration may be necessary.

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1 **MODE**

PUMP	◦ rpm	* Flow Rate
DISP	◦ Time	◦ Volume
PAUSE	◦ Time	◦ PROGRAM
TOTAL	◦	

2 **CAL**

3 **MAX
RESET**

4 **CAL**

1 **MODE**

PUMP	◦ rpm	◦ Flow Rate
DISP	◦ Time	* Volume
PAUSE	◦ Time	◦ PROGRAM
TOTAL	◦	

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Calibrant reset of flow rate

- Change mode to PUMP Flow rate
- Press the CAL button (displayed value blinks)
- Press the MAX/RESET button.
The flow rate is reset to the default value (Default value blinks).
- Confirm with the CAL button
➡ At the same time, the volume is also reset to the default value.

Calibrant reset of volume

- Change mode to DISP Volume
- Press the CAL button (displayed value blinks)
- Press the MAX/RESET button
The volume is reset to the default value (Default value blinks).
- Confirm with the CAL button
➡ At the same time, the flow rate is also reset to the default value.

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1 **MODE**

PUMP	rpm	Flow Rate
DISP	* Time	* Volume
PAUSE	Time	PROGRAM
TOTAL		

▲ ▼

30.0

2 **MODE**

PUMP	rpm	Flow Rate
DISP	* Time	* Volume
PAUSE	Time	PROGRAM
TOTAL		

▲ ▼

800.0

3 **MODE**

PUMP	rpm	Flow Rate
DISP	* Time	* Volume
PAUSE	Time	PROGRAM
TOTAL		

4 **RUN STOP**

a **MAX RESET**

PUMP	rpm	Flow Rate
DISP	* Time	* Volume
PAUSE	Time	PROGRAM
TOTAL		

oder/or/ou

9999 1111

b

PUMP	rpm	Flow Rate
DISP	* Time	* Volume
PAUSE	Time	PROGRAM
TOTAL		

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

Dispensing by volume within a pre-set time

- Change to mode DISP Time. Enter the required time with the ▲▼ keys (0.1 s – 999 h).
 - Change mode to DISP Volume. Enter the required volume in ml with the ▲▼ keys.
 - Push the MODE-button
 - ➔ Both LEDs light up simultaneously.
 - Start pump with RUN/STOP
 - ➔ The pump dispenses the required volume within the pre-set time.
- a. **Changes / Corrections**
Press the Reset key
- ➔ For changing the dispensing time or dispensing volume, MODE returns directly to DISP Time.
 - ➔ Repeat steps 1 to 3
- b. **Error message**
LEDs are blinking alternately
- 9999 = volume too large time too short
1111 = volume too small time too long
- Despite volume and time adjustments, the LEDs keep on blinking:
- ➔ change the pump-head model in order to obtain other flow rate.

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Please note

For entering the number of dispensing cycles see Page 27.

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- ▶ Drehzahl (PUMP rpm)*
- ▶ Fließrate (PUMP Flow rate)*
- ▶ total gefördertes Volumen (TOTAL)

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Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time

1. Go to mode PAUSE Time Enter the pause time (between 0.1 s and 999 h) with ▲▼ keys.

2. Press MODE key once

➔ Both LEDs for PAUSE and DISP Time light up simultaneously. Enter the dispensing time (0.1 s to 999 h) with ▲▼ keys.

3. Start pump with RUN/STOP key

➔ Pump dispenses at intervals until the RUN/STOP button is pushed again.

➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point appears shortly on the display taking effect, however, only with the subsequent dispensing cycle.

➔ During the dispensing process the following values can be retrieved with the MODE-key:

- ▶ rotation speed (PUMP rpm)*
- ▶ flow rate (PUMP Flow rate)*
- ▶ totally dispensed volume (TOTAL)

*The rotation speed and flow rate can be changed with the ▲▼ keys.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time

1. Change mode to PAUSE Time Enter pause time (between 0.1 s and 999 h) with the ▲▼ keys.

2. Press MODE-key twice

⇒ Both LEDs for PAUSE and DISP Volume light up simultaneously Enter the volume with the ▲▼ keys.

3. Start pump with the RUN/STOP key

⇒ Pump dispenses at intervals until the RUN/STOP button is pushed again.

⇒ Active mode blinks.

⇒ With the ▲▼ keys the dispensing volume can be changed even during the dispensing process. The newly entered set point appears shortly on the display taking effect, however, only with the next dispensing cycle.

⇒ During the dispensing process the following values can be retrieved with the MODE-key:

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If necessary, calibrate the volume as indicated on Page 22.

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

- Change into one of the two modes for dispensing at intervals:
 - ▶ DISP Time/Pause, or
 - ▶ DISP Volume/Pause
- Push the CAL-key; the number 0 for an infinite number of dispensing cycles blinks.
 - ➡ With the RESET-key you can switch between 0 and the number of repetitions last entered.
- With the ▲▼ keys the number of dispensing cycles can be adjusted.
- Confirm with the CAL-key.
- Start dispensing at intervals with the RUN/STOP-key.
 - ➡ During the pauses the display shows the remaining number of dispensing cycles.
 - ➡ During long pauses, the remaining interval time can be checked with the CAL-key.
 - ➡ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key; when resuming the process, the pump continues with the remaining amount of dispensing cycles.

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Overcurrent protector

The drive MCP-Z Standard features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL **S** = Speed set-point cannot be reached

OL **I** = Current Overload

Whenever this situation occurs, the pump must be switched off immediately.

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
Pin 2, local / remote
Pin 3, stop / run
Pin 13, speed IN

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Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs

Pin 2, local / remote

For changing between manual control and analog interface. For activating the analog interface, Pin 2 must be connected with Pin 1 (GND).

Pin 3, stop / run

- In normal operation (Pin 2 open) Pin 3 operates the foot-switch (see DIP switch settings, Page 32)
- In remote operation (Pin 2 to GND) the pump starts when connected to Pin 1 (GND)

Pin 5, speed IN

For external speed control (0 – 5 V, 0 – 10 V, 0 – 20 mA, 4 – 20 mA) Alternatives by means of DIP switches inside the pump (see Page 32)

Pin 7, +36 V_{DC}

About +36 V_{DC} are available (max. current 1A).^{*)}

^{*)} The max. total current of the +36 V_{DC} supply on Pin 7 and on Pin 9 and 10 of the valve interface is 1 A. Fuse: F1.6A/250V fast, located on controller board, see diagram on Page 32.

Pin 9, motor speed

The default setting is 0 – 10 V_{DC}, proportionally to the motor speed 60 – 6000 rpm.

Alternatively a frequency range from 0 – 12 kHz is available. Alternatives with sliding switch S2 inside the pump (see Page 32).

Pin 10, +5 V_{DC}

About +5 V_{DC} are available. (max. current 0,5 A)

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Pin 2, local / remote
 Pin 3, stop / run
 Pin 13, speed IN

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Analog interface

Pin 13, speed Pin 5 internal

Depending on how the pump is operated, Pin 13 has different functions:

Analog interface not activated

(Normal operation, i.e. Pin 2 is open) Pin 13 serves as **auto-start function**. If Pin 13 is connected to Pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be ON).

Analog interface activated

(Pin 2 on GND)

- ▶ Pin 13 open: The rotation speed is adjusted via Pin 5 (speed IN).
- ▶ Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
Please note

The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Settings of switch S1

Operation via foot-switch

With DIP switch 5 the user can switch between two possibilities:

- ▶ »FS toggle« (ON/OFF)
- ▶ »FS direct« (ON = as long as pressure is maintained on the footswitch)

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Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5
Pin 3 (footswitch)	FS toggle						OFF*
	FS direct						ON
Pin 5 speed IN	0-5V	470 kΩ		OFF*	OFF*	OFF*	
	0-10V	20 kΩ		OFF	OFF	ON	
	0-20mA	240 Ω		OFF	ON	OFF	
	4-20mA	240 Ω			ON	ON	OFF
Baud Rate RS232	9600 Bd		OFF*				
	1200 Bd		ON				

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Sliding switch S2

This switch affects Pin 9, motor speed

Position A: 0 – 10 V_{DC} (Standard)
Position B: 0 – 12 kHz

➔ **Make sure that the pump is disconnected from the mains supply.**

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➔ Downloads

Valve interface

A detailed description of the valve interface is available for download on our web-site.

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
Valve interface

Up to two valves can be operated directly from the VALVES interface with an adjustable effective voltage between 0 and ca. 26 V, see Page 34, 'valve interface' command section. Commands X,Y,Z,& control the valves. Switching is done with an initial maximum voltage of 26 V_{DC} for 0.5 sec to make sure that the valve moves. The hold voltage can be adjusted with the 'Z_----' and 'X_----' commands through RS232 (but not at the pump's operating panel). A simple terminal program called 'IsmaTerm' is available for PC at www.ismatec.com.


The two 'sense' pins detect if a valve connector is plugged in. If just one valve is used between Pin 1 (valve 2-) and Pin 9 (valve 2+), the 'sense' pins 8 and 14 need to be linked by a piece of wire. If an additional second valve is present at Pin 2 (valve 1-) and Pin 10 (valve 1+), the 'sense' link should be a 47 ohm resistor instead.

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PWM (000 .. 255) connection to GND for valve 1.

Pin 8, sense GND

GND for the detection of a valve plug via short-circuit with Pin 14.

Pin 9, valve 2 +

About +36 V_{DC} are available for valve 2 (max. current 1 A).^{*)}

Pin 10, valve 1+

About +36 V_{DC} are available for valve 1 (max. current 1 A).^{*)}

Pin 14, sense

The valve plug must provide a short-circuit between Pin 8 and 14. This enables the pump to detect the presence of a valve connector.

^{*)} The max. total current of the +36 V_{DC} supply on Pin 7 and on Pin 9 and 10 of the valve interface is 1 A. Fuse: F1.6A/250V, located on controller board, see diagram on Page 32.

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Please note

When switching on the MCP-Z Standard drive, the serial interface will only respond after a delay of 3 seconds.

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
Kommunikationsparameter Communication parameters Paramètres de communication

Standardeinstellung
Default setting
Réglage standard

**9600 baud, 8 bit,
1 stopbit, no parity**
(DIP Switch 1 off)

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Serial interface

RS232 IN (female)

A 9-Pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1 – 8). It is factory-set at 1 and can be changed with @. The address will be stored permanently (even after the pump has been switched off). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed by** a carriage return (ASCII 13). The pump confirms most of the commands with an asterix*. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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Each pump requires an individual address
(see Page 35).

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-Pin RS232 connector.

If the PC is equipped with a 25-Pin connector, you need an additional adapter XC0009.

➔ **Each pump must be allocated its own individual address.**

For this procedure, you first have to connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can you cascade them as shown on the left.

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
Bestell-Nr. XC0009

Programming Software

LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments.

<http://www.ni.com/>

LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our website: <http://www.ismatec.com>

Interface card Moxa C168

including a 25 Pin multiple cable (male) for controlling up to 8 laboratory instruments
Order No. IM0015C168

Connecting cable 9F / 9M

For connecting an MCP-Z Standard to a PC via the 9-Pin RS232 connector. (2 m)
Order No. AG0013

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Pump-head

In general, the MCP-Z Standard gear pump is supplied as a complete unit, consisting of the MCP-Z Standard drive, and as ordered:

1. pump-head
 2. tubing adapters (2 pcs.)
- a. Driving magnet included

For the MCP-Z Standard, a wide range of different gear pump-heads with PTFE, PPS or graphite gears are available.

You find an overview chart of the pump-heads available on Page 49/50.

⚠ Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.

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Mounting the pump-head

Fig. 1

The following tools are required for mounting a pump-head:

- ▶ Allen key 3/32"
- ▶ Phillips screwdriver (Size 0)

Remove all screws from the fixing case (a) of the drive.

(Together with the drive, two long screws are additionally supplied for fastening the large pump-heads of the series 200 (at the top and the bottom) of the pump-head face).

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Fig. 1

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Mounting the pump-head

Fig. 2

Fasten the driving magnet on the motor shaft (through the hole underneath).

➔ Ensure that there is a gap of about 3 mm between the magnet and the brim of the fixing case.

⚠ **During operation the back of the magnet must not touch the housing of the pump drive.**

Fig. 3

Screw the two tubing adapters into the pump-head.

➔ Before screwing the tubing adapters into the pump-head, we recommend you wrap **PTFE tape around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.

Fig. 4

Screw the pump-head to the fixing-case of the drive. The face of the pump-head with the engraved model number must be directed to the user. The engraved arrow indicates the flow direction.

➔ Please note that gear pumps can only pump in one direction. An arrow engraved on the front of the pump-head indicates the correct flow direction.

⚠ **Fill the pump-head with liquid before you start pumping. Avoid dry-running of the pump-head.**

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Fig. 2

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- ▶ PTFE = Polytetrafluorethylen
- ▶ PPS = Polyphenylensulfid
- ▶ Graphite = Kohlenstoff
- ▶ PEEK = Polyätherätherketon

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Gear Pump Technology

»Cavity style« pump-head

- ▶ The »Cavity style« pump-head can be used for viscous media and for applications with a certain suction lift
- ▶ Excellent chemical resistance
- ▶ Smooth operation at a low noise
- ▶ Low inner friction
- ▶ Incl. B mount (flange)

»Suction shoe« pump-heads

- ▶ Enhanced pump performance at elevated differential pressures
- ▶ Fewer parts, easy to service
- ▶ Suitable for elevated temperature range
- ▶ Not suitable for applications requiring a suction lift
- ▶ Incl. B mount (flange)

Pump-head specifications

- ▶ Stainless steel pump-head housing SS 316/DIN 1.4571

Temperature range

- ▶ -46°C to +54°C / -51°F to +129°F (Cavity style)
 - ▶ -46°C to +177°C / -51°F to +350°F (Suction shoe)
- (higher temperatures possible on request)

Differential pressure

The difference between the pressure on the inlet and the discharge port of the pump-head max. 5.2 bar (75 psi) depending on pump-head (Page 49/50)

System pressure

The pressure within the pump-head max. 21 bar / 305 psi (Exception: Z-1830, Z-150: 22 bar)

Gear material

- ▶ PTFE = Polytetrafluoroethylene
- ▶ PPS = Polyphenylenesulphide
- ▶ Graphite = Carbon graphite
- ▶ PEEK = Polyetheretherketone

Ports

(Inside thread) 1/8"-27 NPT for tubing with 6.4 mm ID

The pump-heads are interchangeable. Different drives take different magnets. The MCP-Z Standard drive is originally supplied with the correct magnet.

When ordering a new pump-head, please state the intended operation temperature. If you intend to use a Cavity style pump-head for temperatures exceeding 54°C (130°F), the correct sealing material must be ordered. Please contact your nearest ISMATEC® representative.

Some pump-heads are available with an **internal bypass**. This bypass is to be used for safety reasons only and must not be used as pressure control (see Page 28).

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Pumpenkopf Nr. Pump-head No. Tête de pompe No.	181	183	183HC	186	186P	186HC	1830	1830 P	200	200-P	200 PKC	201	201 PKC
Bestell-Nr. Order No. No. de commande	MI 0007	MI 0008	MI 0310	MI 0006	MI 0312	MI 0309	MI 0131	MI 0280	MI 0022	MI 0306	MI 0341	MI 0023	MI 0378
Fließrate ml/min max. ¹⁾ Flow rate Débit	252	504	540	102	102	102	550	550	3380	3380	3380	7020	7020
Fließrate ml/min max. ¹⁾ Flow rate Débit	2.5	5.0	5.0	1.0	1.0	1.0	5.5	5.5	34	34	34	70	70
Differenzdruck max. (bar) Pression différentielle (psi) Differential pressure	2.8 40.6	2.8 40.6	2.8 40.6	1.4 20.3	2.3 33.3	1.4 20.3	5.2 75.4	5.2 75.4	5.2 75.4	5.2 75.4	5.2 75.4	3.5 50.8	3.5 50.8
Zahnradmaterial Gear material Matériaux engrenage	Graphite	Graphite	Graphite	Graphite	PEEK	Graphite	PPS	PEEK	PPS	PEEK	PEEK	PPS	PEEK
Dichtungen Seals Etanchéités	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	Viton®	Viton®	Kalrez® **)	Viton®	Kalrez® **)
Gehäusematerial Housing material Matériaux caisse	SS316	SS316	*) Hastelloy C276	SS316	SS316	*) Hastelloy C276	SS316	SS316	SS316	SS316	SS316	SS316	SS316
Systemdruck max. (bar) System pressure Pression de système (psi)	21 305	21 305	21 305	21 305	21 305	21 305	21 305	21 305	21 305	21 305	21 305	21 305	21 305
Betriebstemperatur °C Operating temperature Temperature operation °F	-46...+177 -50...+350	-46...+177- -50...+350	-46...+177 -50...+350	-46...+177 -50...+350	-466...+177 -50...+350	-46...+177 -50...+350	-46...+177 -50...+350	-46...+177 -50...+350	-29...+177 -20.2...+350	-29...+177 -20.2...+350	-29...+177 -20.2...+350	-29...+177 -20.2...+350	-29...+177 -20.2...+350
Bypass Internal bypass Bypass interne	--	--	--	--	--	--	--	--	✓	--	✓	✓	✓
Anschlüsse (Innengewinde) Ports (internal thread) Ports (internal thread)	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT

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Pumpenkopf Nr. Pump-head No. Tête de pompe No.	120	130	Z-140	140 P	140 HC	142	142 HC	150	150 WI
Bestell-Nr. Order No. No. de commande	MI 0013	MI 0019	MI 0016	MI 0313	MI 0284	MI 0018	MI 0311	MI 0020	MI 0265
Fließrate ml/min max. ¹⁾ Flow rate Débit	3840	3840	3840	3840	3840	5640	5640	3840	3840
Fließrate ml/min max. ¹⁾ Flow rate Débit	38	38	38	38	38	56	56	38	38
Differenzdruck max. (bar) Pression différentielle (psi) Differential pressure	3.5 50.8	5.2 75.4	3.5 50.8	5.6 81.2	3.5 50.8	3.5 50.8	3.5 50.8	5.2 75.4	5.2 75.4
Zahnradmaterial Gear material Matériaux engrenage	PTFE	PPS	PTFE	PEEK	PTFE	PTFE	PTFE	PPS	Nickel-Carbid
Dichtungen Seals Étanchéités	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE
Gehäusematerial Housing material Matériaux caisse	SS316	SS316	SS316	SS316	** ²⁾ Hastelloy C276	SS316	** ²⁾ Hastelloy C276	SS316	*** ³⁾ Oberfläche Surface hardened
Systemdruck max. (bar) System pressure Pression de système (psi)	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6
Betriebstemperatur °C Operating temperature Temperature operation °F	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129
Bypass Internal bypass Bypass interne	✓	✓	--	--	--	--	--	--	--
Anschlüsse (Innengewinde) Ports (internal thread) Ports (internal thread)	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT

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Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

We strongly recommend you thoroughly flush the pump-head with water after use.

Provided the MCP-Z Standard is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear by himself (see Page 45).

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Repairs

For repairs please send the MCP-Z Standard drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., model.

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Spare parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- ▶ spare parts
- ▶ parts lists
- ▶ wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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Pump is not running.

Is the pump connected to the mains?

Is the voltage selector set correctly? Are the correct fuses inserted?

➔ see Pages 8 – 9

Flow rate is incorrect.

Is the correct pump-head ID entered in the actual program?

Have you calibrated the flow rate or dispensing volume in the actual program?

➔ see Pages 16/19/22

No response from the RS232 interface?

Is the pump address correct? Is the pump connected to the correct PC port (COM1,...)?

Do the transmission parameter of the pump and the PC correspond with each other (Baud rate, ...)?

➔ see Pages 33 – 41

Cascaded pumps do not respond.

Are the transmission parameters of the pump and PC correct (Baud rate, ...)?

Have all the connected pumps the correct address?

➔ see Page 42

Should problems arise that are not listed here, please contact your local ISMATEC® distributor.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Drive

Motor type	DC-Motor
Speed range	60 – 6000 rpm
digitally adjustable	in steps of 1 rpm

Differential pressure

max. 5.6 bar (81.2 psi) depending on pump-head see Pages 49 – 50

Remote control

via RS232 and analog interface (see Pages 30 – 41)

Mains connection

230 V_{AC} 50 Hz
115 V_{AC} 60 Hz

Fuse rating

(230 V) 2 x T1.25A/250V
(115 V) 2 x T2.5A/250V

Power consumption

Max. 150 W

Operating conditions

Temperature	41 to 104°F,
Rel. humidity	max. 80%
– not condensing, at normal laboratory conditions	

Dimensions/Weight (Drive)

DxWxH	8.7 x 6.1 x 10.3 inch
Weight (without pump-head)	14.1 lb

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


Centering flange

The centering flange fits onto MCP-Z drive. It is only required for the Micropump 200 and 201 pump-head series where it facilitates their precise mounting.

Accessories:

- a. 4 flat-headed screws
M 2.5 x 8 BN0187
- b. 1 centering
flange IS10379

Mounting instructions:

1. Please ensure that the circular drive magnet (c) is correctly mounted.
2. When mounting the centering flange, make sure that
 - ▶ its elevated inner ring (1) exactly fits into the drive flange, and
 - ▶ the two straight holes (2) are placed in a vertical position at the top and the bottom.

➡ When screwing the centering flange onto the drive, make certain that none of the four flat-headed screws protrude.
3. You can now proceed with mounting the gear pump-head with the two screws in the usual way.

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ISMATEC®

Learn more
about the
Ismatec® product
line by visiting:
www.ismatec.com

For ordering and technical support, please contact:

North America sales@ismatec.com | 1-800-323-4340 | 1-847-549-7600

Europe sales.europe@ismatec.com | +49 (0) 9377 9203-0

 **Cole-Parmer®**