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FMI007 Keramikkolben

mit Settings-Menü

Valveless Rotary Piston
Dispensing Pump
Microprocessor controlled

Drive
ISM321

Head 2–25 µl/stroke

FMI009 ceramic piston
FMI010 steel piston

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Piston pump-heads require
thorough flushing after use

It is recommended that the pump-head be flushed at a small back pressure **after every pumping process** to prevent the piston from seizing due to remaining particles (especially after pumping media containing salt, protein or particles).

Do not pump fluids containing suspended solids that are abrasive.

Please observe our note regarding filtration in the chapters Safety Precautions and Maintenance.

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Caution

Rotating pumphead creates a pinch and crush hazard.

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Safety precautions

ISMATEC® ceramic piston pumps are designed for pumping and dispensing applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- ▶ The pump must not be operated outside the designed operating and environmental conditions. In particular, it is the user's responsibility not to exceed the specified maximum differential pressure for the pump-head.
- ▶ The pump must not be used:
 - for medical applications on human beings,
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- ▶ The circuit between mains supply and pump has to be connected to earth ground.
- ▶ Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- ▶ The pump must be switched off before taking off the pump-head.
- ▶ Do not open or remove the housing while the pump is operating.

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Media containing particles:

It is recommended that the pump head be rinsed **after every pumping process** to prevent the piston from seizing due to remaining particles (especially after pumping media containing salt, protein or particles).

Do not pump media containing particles that exceed a diameter of 0.8 mm.

Piston stroke or stroke volume must be set large enough with respect to the particle size!

Piston fractures are not covered by our warranty.

Corrosive liquids:

Corrosive liquids may eventually produce leak paths around the sealing surfaces of the pump-head. The pump should be inspected for leaks on a regular basis.

Dry running:

Before starting to run the pump, it is advisable to fill the pump-head with liquid.

Ensure that while priming, the pump does not run dry for more than 15 seconds.

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Safety precautions

- ▶ Only new fuses, according to the specifications stated on Page 10 in this manual, must be used.
- ▶ The fuse-holder must not be short-circuited.
- ▶ The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when routing tubing in explosion-proof chambers.
- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- ▶ Repairs may only be carried out by a skilled person who is aware of the potential hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Pump drive: 2 years

Pump-head

We advise you read the enclosed operating manual of the pump-head manufacturer.

Piston fractures are not covered by our warranty (see our instructions for flushing the pump-head on Page 4)

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions. If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packaging or in a packaging of equal quality. For pump-heads from manufacturers other than ISMATEC® the warranty terms of the specific manufacturer are valid. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use
- ▶ unauthorized modification or misuse by the user or by a third-party
- ▶ improper site preparation and maintenance
- ▶ operation outside of the environmental and electrical specifications for the product
- ▶ use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

The REGLO-CPF *Digital* is a valveless rotary piston dispensing drive.

Package:

- ▶ Pump drive ISM321
- ▶ Power Cord
- ▶ Operating Manual
- ▶ Pump-head as ordered

Head 2–25 µl/stroke

FMI009	ceramic piston
FMI010	steel piston

Head 5–50 µl/stroke

FMI005	ceramic piston
--------	----------------

Head 10–100 µl/stroke

FMI007	ceramic piston
--------	----------------

Please check the package

and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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Rear panel

1. RS232 IN (female) for operation controlling via PC
2. RS232 OUT (male) for cascading up to 8 pumps

Digital input (TTL level)
Run/Stop, autostart

Digital output 0–10 kHz (speed output)

3. Mains socket
4. Fuse-holder

⚠ Mains voltage

Mains voltage	Fuse rating
100–240 V _{AC}	1 x T2A/250V
50/60 Hz	

⚠ Socket/Power Cord

Use exclusively the originally supplied power cord.

The socket must be connected to earth ground (protective conductor contact).

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Changing the fuses



➡ Switch the pump OFF, pull out the mains plug.


1. Pull out the fuse-holder by opening it with a small screw-driver (size 0).
2. Remove fuse (a) and replace it with the spare fuse (b) 1 x T2A/250V.

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Operating panel (* = Settings Page 13)

1. ON/OFF switch
2. Digital LED display
3. Increase **value**
4. Reduce **value**
*settings
5. **RUN/STOP** (double function)
 - a. Starts and stops the pump
 - b. Interrupts and continues a dispensing cycle (Page 23)

➔ The function Start/Stop can also be triggered by an external device such as a footswitch (accessories on Page 43).

6. Double function

a. Rotation direction

Changes the rotation direction. By changing the rotation direction or starting the pump, the rotation direction is indicated in the display with a running hyphen.

b. RESET

- Resets the accumulated volume in mode TOTAL to zero (Page 20)
- Interrupts dispensing cycle for good (Page 23)

*escape

7. Double function

- a. **CAL** = Calibrating key for the flow rate or dispensing volume
- b. **MAX** = Max. speed in the modes PUMP•Flow rate and TOTAL (keep the MAX key pressed down while the pump is running)
*ok

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

Operating panel (* = Settings Page 13)

8. MODE key Selecting operating mode
* default

9. LEDs for active operating mode

- ▶ **PUMP • Flow rate**
Pumping by flow rate (ml/min)
- ▶ **DISP • Volume**
Dispensing by volume (ml)
- ▶ **DISP • Time**
Dispensing by time
- ▶ **PAUSE • Time**
Pause time for dispensing at intervals
- ▶ **TOTAL**
Puming with read-out of totally
delivered volume

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Start-up information

After switching the power supply switch
ON, the following values are displayed:

1. LED test »8.8.8.8.«
2. Firmware version
3. Displays the operating mode last used

⚠ **Before starting the pump for
the first time or after connecting
a new pump-head,** the ID-code of
the pump-head must be entered
(see basic settings, Page 17).

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Please note

The settings menu requires a break of 5 seconds between power off and on.

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Start-up information


1. Switch the pump OFF and wait for 5 seconds.
2. Keep the settings key pressed and switch the pump ON.
3. Select the basic settings by using the ▲▼ keys
 - Cycles (dispensing cycles)
 - Back-steps (number of back-steps for volume dispensing)
 - Address (RS232, pump address)
 - Footswitch (Toggle/Direct)
 - Autostart (Standby/RUN)
 - Digits (Number of digits displayed for flow rate and dispensing volume)
 - Display Speed (Type of speed display)
 - Stroke volume
 - Head = ID-Code (pump-head)
 - Init (resets all basic settings to the default values)
4. **ok** = confirm your selection
5. Select the value by using the ▲▼ keys, or
 - press escape for cancelling wrong input
 - the values can be reset to the default value by using the default key
6. **ok** = confirm your selection
7. escape
Quit menu basic settings

➔ Basic settings refer to the following Pages

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Number of dispensing steps 0..9999
(0 = Default value, infinity)

Backsteps

Number of stroke volume back-steps (for drip-free dispensing in the mode DISP Volume) from 0 .. 100 (0 = Default value)

Address (RS232)

Pump address 1-8 (Page 34)
(1=default value)

Footswitch

- ▶ Toggle (default value)
 - Step on footswitch = pump starts running
 - Press footswitch again = pump stops
- ▶ Direct
the pump runs as long as the footswitch is pressed.

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Head Identification code (ID) of the pump-head, see also table on Page 42 (25 = default value)

➔ In order to obtain correct pumping and dispensing values, the individual ID code of the mounted pump-head must be entered; especially when a pump is used for the first time and each time the pump-head is changed.

Stroke volume

Configuring the stroke volume. Choose a value between 1µl and the max. stroke volume of the mounted pump-head (25, 50 or 100 µl).

Confirm the adjusted stroke volume with the **ok** key.

➔ Turn the knurled adjustment ring until the value on the fixed and revolving flow control ring equals the digital LED read-out. (The fixing screw must not be loosened!)

➔ The text SET/HEAD/TO is repeated five times. (skip repetitions with CAL/MAX-key).

➔ For accurate dispensing, it is recommended that the stroke volume not be set below 10% of the max. value.

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
Basic settings

Init (Initializing)

Pressing the **ok** key resets the following basic settings to the default values:

Mode	PUMP•Flow rate
Speed	400 rpm
Dispensing volume	54 revs
Dispensing time	4.50 sec
Pause time	2.00 sec
Rotation direction	clockwise
Back-steps	0
Number of dispensing cycles	0
Pump-head digits display	3
Footswitch	Toggle
Autostart	Standby
Display	Flow rate

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Zero point calibration

Before using the pump the first time, or whenever necessary, it is recommended that the 0-point of the flow control ring be calibrated.

- ▶ Start the pump in the mode PUMP • Flow rate with the RUN/STOP key and pump the liquid to be delivered.
- ▶ While the pump is in operation, turn the knurled adjustment ring to the right until the pump stops delivering liquid (the motor is still running!).
- ▶ Loosen the fixing screw of the revolving flow control ring and turn it to the value '0'.
- ▶ Tighten the fixing screw again.
- ▶ Adjust the required stroke volume and configure the basic setting stroke volume accordingly (see Page 17).

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Enter the stroke volume in the settings menu (see basic settings, Page 17).

The flow rate is determined by the stroke volume and the rotation speed.

➔ Flow rate = Stroke volume x speed

However, there is no recommendation of how to set the stroke volume and rotation speed in an optimum way.

Basically the following applies:

- ▶ *Small stroke volume/ high speed:*
 - Slow pulsation
 - better resolution of the flow rate and dispensing volume
 - tendency to cavitation (insufficient flow-through of the liquid in the pump-head)
- ▶ *Large stroke volume/ low speed:*
 - for elevated viscosity
 - for differential pressures
 - for reduced resolution of the flow rate and dispensing volume

Important: The pump always dispenses integer multiples of the stroke volume. If the stroke volume is adjusted to 10 µl and the targeted volume is 75 µl, the pump permits only 70 or 80 µl when selecting the dispensing volume with the ▲▼ keys.

➔ Adjust stroke volume mechanically to 7.5 µl, i.e. for a 25 µl pump-head from 180 units to 135 units.

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Pumping by flow rate

Before entering the required flow rate, please ensure that the “Display Speed” in basic settings is set to “Flow” (flow rate in ml/min) (see Pages 16 and 14). If set to “rpm” (rotation speed), neither can the flow rate be entered nor a calibration be carried out.

- Change mode to PUMP • Flow rate.
- Enter the required flow rate with the ▲▼ keys (displayed in ml/min)

Start pump with RUN/STOP

➔ The flow rate can also be adjusted while the pump is running.

➔ For an accurate flow rate, it is recommended that the pump be calibrated (see Page 21).

- For fast filling or emptying the system, keep pressing the CAL/MAX key while the pump is running.

➔ During the pumping process the following values can be retrieved with the MODE key:

- Flow rate (PUMP rpm)
- Totally pumped volume (TOTAL)

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TOTAL volume

Change with the MODE key to TOTAL. The totally delivered volume is displayed.

Four display readings are available:

Microlitre:	e.g.	17.0 µl	=	17.0µ
Millilitre:	e.g.	17.0 ml	=	17.0
Litre:	e.g.	17.0 l	=	17.0L
Litre:	e.g.	1700 l	=	1700.

For zero-setting, either

- press the RESET key once, or
- switch the pump OFF.

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Before performing a calibration, please ensure that the "Display Speed" in basic settings is set to "Flow" (flow rate in ml/min) (see Page 16 and 14). If set to "rpm" (rotation speed), neither can the flow rate be entered nor a calibration be carried out.

- Change mode to PUMP • Flow rate.

Set the required flow rate with the ▲▼ keys

➔ Check: have you configured the required stroke volume and adjusted the revolving flow control ring accordingly? (see Page 17)

- Change the mode to DISP•Time

– Enter 60 secs. with the ▲▼ keys
 – Start the pump with the RUN/STOP key

➔ The pump stops automatically after 60 secs.

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Calibrate the ascertained value as follows:

- Return to the mode PUMP • Flow rate

- Press the CAL/MAX key (displayed value blinks)

Enter the measured value (balance) with the ▲▼ keys

- Save with the CAL/MAX key. The pump displays the value to be adjusted with the flow control rings (0,0–450,0 div.) The text SET/HEAD/TO is repeated five times. (Skip repetitions with CAL/MAX-key).

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Calibrating the flow rate

➔ Turn the knurled adjustment ring until the correct value on the fixed (red) and revolving flow control ring (black) equals the digital LED read-out. A complete revolution equals 100 divisions (max. 450 divisions).

The fixing screw must not be loosened!

6. The pump is now calibrated and ready. Start the pump with the RUN/STOP-key.

➔ In case of problems set the default calibration first (see Page 25).

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Please Note

A dispensing cycle can be interrupted by pressing the RUN/STOP key.

- the volume dispensed up to now is blinking
- by pressing the RUN/STOP key again, the dispensing cycle is continued
- by pressing the RESET-key the dispensing cycle is stopped for good

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Dispensing by volume

1. Change mode to DISP • Volume
2. Use the ▲▼ keys for entering the required dispensing volume. Two display readings are available:
 Microlitre: e.g. 12.5 µl = **12.5µ**
 Millilitre: e.g. 230.5 ml = **230.5**
 ➔ For dispensing accurately, it is recommended that the pump be calibrated (see Page 24).
3. Start the pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

- ➔ The dispensing speed can be adjusted in the mode PUMP • Flow rate.
- ➔ With the ▲▼ keys the dispensing volume can also be changed even during the dispensing process. The newly entered set point appears briefly in the display.
- ➔ During the dispensing process the flow rate (PUMP • Flow rate) and the total volume can be retrieved by pressing the MODE key.
- ➔ The flow rate can be changed even during the dispensing process by using the ▲▼ keys.

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Calibrating the volume

1. Change the MODE to DISP • Volume

➔ Use the ▲▼ keys for entering the required dispensing volume.

2. Start with RUN/STOP

➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes.

➔ Measure the dispensed liquid by volume or weight.

3. Stay in the same mode (DISP • Volume) and press the CAL/MAX-key (display blinks)

4. Enter the effectively dispensed volume with the ▲▼ keys.

5. Save with the CAL/MAX key. The pump displays the target value for the adjustment. The text SET/HEAD/TO is repeated five times.

➔ Turn the knurled adjustment ring till the correct value on the fixed (red) and revolving flow control ring (black) equals the digital LED readout. A complete revolution equals 100 divisions (max. 450 divisions). The fixing screw must not be loosened!

➔ The pump is now ready for calibrated dispensing.

6. Start the pump with RUN/STOP

➔ Depending on the application conditions repeated calibration may be necessary.

➔ A dosing volume of less than 10 strokes can not be calibrated.

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Default calibration of flow rate

1. Change mode to PUMP • Flow rate.
2. Press the CAL/MAX key (displayed value blinks)
3. Press the **default** key. (Default value blinks)
4. Confirm with the CAL/MAX key. The pump displays the value to be adjusted with the flow control rings (0.0–450.0 div.). The fixing screw must not be loosened!

➔ At the same time, the volume is also reset to the default value.

The text SET/HEAD/TO is repeated five times. (Skip repetitions with CAL/MAX key.)

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Default calibration of volume

1. Change mode to DISP • Volume
2. Press the CAL/MAX key (displayed value blinks)
3. Press the default key (Default value blinks)
4. Confirm with the CAL/MAX key. The pump displays the value to be adjusted with the flow control rings (0.0–450.0 div.). The fixing screw must not be loosened!

➔ At the same time, the flow rate is also reset to the default value.

The text SET/HEAD/TO is repeated five times. (skip repetitions with CAL/MAX-key).

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Dispensing by volume

The dispensing time can be entered from 0.1s to 999h.

1. Change the mode to DISP•Time
2. Enter the required dispensing time with the ▲▼ keys
 ➔ The display accelerates when pressure on the ▲ or ▼ maintained
 Seconds: 0.1–899.9 (in 0.1s steps)
 Minutes: 15'–899' (in 1min steps)
 Hours: 15h–999h (in 1h steps)
3. Start the pump with RUN/STOP

➔ With the ▲▼ keys the dispensing time can be changed while the dispensing cycle is in process.

➔ During the dispensing process the flow rate (PUMP • Flow rate) and the total volume can be retrieved by pressing the MODE-key. The flow rate can be changed during the dispensing process by using the ▲▼ keys.

➔ Very short dispensing times can result in non-repeatable dispensing volumes.

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



Dispensing by volume within a pre-set time

1. Change to mode DISP • Volume. Enter the required volume (in µl or ml) with the ▲▼ keys.
 2. Change to the mode DISP•Time. Enter the required time with the ▲▼ keys (0.1s – 999h, see Page 26)
 3. Push the MODE-key once
 ➔ Both LEDs light up simultaneously
 ➔ 1111 and 9999 indicate an error message (see below).
 4. x Start pump with RUN/STOP
 ➔ The pump dispenses the required volume within the pre-set time.
- a. Changes / Corrections**
 For changing the dispensing time or dispensing volume press the RESET-key when the pump has stopped.
 ➔ The pump returns directly to MODE DISP • Volume
 ➔ Repeat steps 1 to 3
- b. Error message**
 LEDs are blinking alternately
- 9999 = volume too large
 time too short**
- 1111 = volume too small
 time too long**

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Please note

For entering the number of dispensing cycles see Page 15.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time.

- Go to the mode PAUSE • Time. Enter the pause time (between 0.1s and 999h) with ▲▼ keys.
- Press MODE key until both LEDs for PAUSE and DISP • Volume light up simultaneously.
➔ Enter the dispensing volume with the ▲▼ keys.
- Start with the RUN/STOP key.
➔ The pump dispenses now at intervals.
➔ **Active** mode blinks.
- Pushing the RUN/STOP key again interrupts the dispensing process or pause period.
➔ The display blinks.

The user can now proceed as follows:

- ➔ By pressing the RUN/STOP key again, the dispensing or pause period is continued.
- ➔ By pressing the RESET-key the pump is stopped for good.
- ➔ With the ▲▼ keys the dispensing volume or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.
- ➔ During the dispensing process the flow rate can be retrieved with the MODE key and adjusted with the ▲▼ keys. Furthermore, the total volume is displayed.

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Please note
For entering the number of dispensing cycles see Page 15.

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**Intermittent dispensing
(by time)**

Intermittent dispensing by time with a pre-set pause time.

- Go to mode PAUSE • Time Enter the pause time (between 0.1s and 999h) with ▲▼ keys.
 - Press MODE key until both LEDs for PAUSE and DISP•Time light up simultaneously.
➔ Enter the dispensing time (0.1s to 999h) with ▲▼ keys.
 - Start with the RUN/STOP key.
➔ The pump dispenses now at intervals.
➔ Active mode blinks.
 - Pushing the RUN/STOP key again interrupts the dispensing process or pause period.
➔ The display blinks.
- The user can now proceed as follows:
➔ By pressing the RUN/STOP key again, the dispensing or pause period is continued.
➔ By pressing the RESET-key the pump is stopped for good.
- ➔ With the ▲▼ keys the dispensing or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.
- ➔ During the dispensing process the flow rate can be retrieved with the MODE key and adjusted with the ▲▼ keys.

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Please note

With RUN/STOP a dispensing cycle can be interrupted and resumed.

With RESET an interrupted dispensing cycle can be definitely stopped.

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

See basic settings on Page 15.

➔ During the pauses the display shows the remaining number of dispensing cycles.

➔ During long pauses, the remaining interval time can be checked with the CAL/MAX key and adjusted by using the ▲▼ keys.

➔ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key. When resuming the process, the pump continues with the remaining number of dispensing cycles.

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Drip-free dispensing

Reverse piston strokes allow drip-free dispensing (1 to 100 strokes).

See basic settings bStP (backsteps) Page 15.

Reverse piston strokes are active in the following modes:

- Dispensing by volume (see Page 23)
- Dispensing a volume within a pre-set time (see Page 27)
- Intermittent dispensing by volume (see Page 28)

 **Reject the first dispensing step. Afterwards, the number of reverse piston strokes are compensated.**

Deleting the reverse piston strokes
See basic settings CYCL, Page 15.

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Pumping against pressurerate

The REGLO-CPF *Digital* can be used for continuous duty at a differential pressure of max. 6.9 bar.

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When the pump is idle

It is recommended that the pump head be flushed **after every pumping process** to prevent the piston from seizing due to remaining particles (especially after pumping media containing salt, protein or particles).

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Overcurrent protector

The drive REGLO-CPF *Digital* features an overload protector. It disables the pump drive in case of a motor blockage. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

Whenever this situation occurs, **the pump must be switched OFF immediately**. Wait for 5 sec. before switching ON the pump again. Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure etc.).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
Please note

The analog interface is integrated in the port 'RS232 OUT' and can be used for remote control (e.g. by footswitch).

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Analog interface

Pin 4: +5VDC

+5VDC For connecting an external control unit (current <100mA).

Pin 5: GND

GND Point of reference for all signals described.

Pin 6: start

start Connection for an external trigger device (footswitch, etc.) For operating modes, refer to Basic settings "Footswitch" Page 15.

Pin 8: busy

busy Informs about pump status.
Pumpe in operation
= output +5V
Pumpe in stand-by
= output on ground
(The output resistance is 4.7 k Ω)

Pin 9: speed OUT

speed OUT A square-wave signal with a frequency proportional to the speed (rpm) is at your disposal:
120 pulses are generated per one rotor revolution.

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
Please note

When switching on the REGLO-CPF Digital pump, the serial interface will respond after a delay of approximately 5 seconds.

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Serial interface

RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1–8). It is factory-set at 1 and can be changed with @. The address will be stored permanently (even after the pump has been switched off). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string is **completed** by a carriage return (ASCII 13). The pump confirms most of the commands with an asterisk *. Yes/No inquiries are answered by + (yes) or – (no). Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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
RS232 OUT (male)

Pin 2: RS232 Rx

Is used for connecting additional pumps together with Pin 3/5.

Pin 3: RS232 Tx

Is used for connecting additional pumps together with Pin 2/5.

Pin 5: GND

Point of reference for all signals described.

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Each pump requires
an individual address.
(see Page 34)

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009 .

➔ **Each pump must be allocated its own individual address.**

For this procedure, first connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can they be cascaded as shown on the left.

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Please note

In order to reach the max. flow rate, the stroke volume must be set to the max. value (25, 50 or 100 µl) depending on the pump-head model.

If, depending on the application, this stroke volume must be reduced (see Pages 17–18), it is recommended that settings not be below 10% of the maximum value.

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Footswitch

Order No. ISM894

This footswitch serves as a Start/Stop device. It is very useful when using the pump as a dispenser for filing tubes, bottles, etc. Both hands are free for handling bottles etc.

Connecting cable 9F / 9M

For connecting a REGLO-Z Digital to a PC via the 9-pin RS232 connector. (2 m)
Order No. AG0013

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Maintenance

It is recommended that the pump-head be flushed thoroughly after pumping in order to prevent the piston from seizing due to remaining particles (especially when pumping media containing salt, proteins or particles).

Otherwise, the REGLO-CPF Digital requires no particular maintenance or mechanical adjustments.

Avoid pumping media containing abrasive solids.

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Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Please note

For maintenance, service, repairs and warranty of the pump-head, refer to the enclosed instruction manual of the pump-head manufacturer. In case of doubt, that manual shall prevail.

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Service and repairs

Like all pumps, the REGLO-CPF *Digital* contains components that will wear over a period of time.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packaging or a packaging of equal quality.

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Replacement parts

For repairs performed by the owner (out of the warranty period), ask your ISMATEC® distributor for:

► Piston-/cylinder group

Order no.	Type
FMI311	H-00CVC-LF0.50
FMI310	H-00SKY-LF
FMI307	H-0CKC
FMI308	H-1CKC

► Lip seals (2 Prices are needed)

Order no.	Type
FMI325	H408-00J (RH00)
FMI312	H408-0A (RH0)
FMI319	H408-1A (RH1)

- Parts lists
- Replacement parts for pump drive
- Wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Speed
40–1800 rpm, digitally adjustable in steps of 1 rpm

Differential pressure
max. 6.9 bar (100 psi)

Remote control
via RS232 and analog interface (see Pages 33–41)

Mains connection
100-240 V_{AC}/50/60Hz

Fuse rating power supply
1x T2A/250V

Power consumption
Max. 140 W

Operating conditions
Temperature +41 to +104°F
Rel. humidity max. 80%
– not condensing, at normal laboratory conditions

Dimensions/Weight
DxWxH 9.8x3.9x5.3 inch
Weight 4.6 lb

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
Assembly parts from ISMATEC®

1. REGLO-CPF drive
2. Coupling completely assembled
3. Mounting plate
4. Phillips Flat-head Machine Screws (4 pcs.)
5. Allen key 2.5 mm from FMI

from FMI


- a. CPF pump-head (RH00, RH0, RH1)
- b. Screws (2 pieces)
- c. Serrated lock washers (2 pieces)

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
- 1. Attaching the mounting plate**
Remove screws (b) and serrated lock washers (c) from the pump-head.
- Fasten the mounting-plate (3) on the pump-head with the 2 screws (b) and serrated lock washers (c).

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
- The 4 countersinks (5) must be concave to the pump-head.
- 3. Fixing the coupling**
Place the completely assembled coupling (2) onto the pump-head shaft and push it on as far as it will go. The Allan screw (d) must be centered on the flat of the pump-head shaft (e).
 - Tighten the coupling (2) with the Allan key size 2.5 mm (5).

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


- ▶ The scale of the red flow control ring (f) must be pointing upwards
- ▶ Hold the mounting plate in such a way, that the screws (b) are on the top and at the bottom
- ▶ Turn the coupling, so that afterwards the Allan screw (d) can be tightened through the aperture (g) with the Allan key (5)
- ▶ Insert the pump-head assembly and push it on as far as it will go. Check the position of the Allan screw at the aperture (g).

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Mounting a ceramic piston pump-head on the REGLO-CPF drive

- Drive in the Flat-head Machine Screws (4). Tighten the Allan screw (d) with an Allan key (5) size 2.5 mm.
- Configuration of the pump-head in the basic settings (Page 17).

Head Identification code (ID) of the pump-head, see also table on Page 42. (25 = default value).

➔ In order to obtain correct values, the individual ID code of the mounted pump-head must be entered; especially when a pump is used for the first time and each time the pump-head is changed.

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ISMATEC®

Learn more
about the
Ismatec® product
line by visiting:
www.ismatec.com

For ordering and technical support, please contact:

North America sales@ismatec.com | 1-800-323-4340 | 1-847-549-7600

Europe sales.europe@ismatec.com | +49 (0) 9377 9203-0

 **Cole-Parmer®**