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## Gear pump with variable speed drive

### **REGLO-Z Analog**

Drive without pump-head  
Order No. ISM895

### **REGLO-ZS Analog**

Control unit with separate  
drive, without pump-head  
Order No. ISM896

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### Caution

Rotating pumphead creates a pinch and crush hazard.

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## Contents

|                               |       |
|-------------------------------|-------|
| Safety precautions            | 3–5   |
| Warranty terms                | 6–7   |
| Product                       | 8     |
| Operating panel               | 9     |
| Rear panel                    | 10    |
| Mains voltage                 | 10–11 |
| Voltage setting               | 11    |
| Changing the fuses            | 11    |
| Starting the pump             | 12–13 |
| Analog interface              | 14–16 |
| Pumping against pressure      | 17    |
| Before switching the pump off | 17    |
| Pump-head                     | 18    |
| Mounting the pump-head        | 19–20 |
| Gear Pump Technology          | 21    |
| Pump-head Specifications      | 22–25 |
| Pump-head with bypass         | 26    |
| Accessories                   | 27–28 |
| Maintenance, Repairs          | 29    |
| Spare parts, Disposal         | 30    |
| Technical specifications      | 31    |

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## Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- ▶ The circuit between mains supply and pump has to be connected to earth ground.
- ▶ The pump must not be operated outside the designed operating and environmental conditions.
- ▶ The pump **must not** be used:
  - for medical applications on human beings
  - in explosion proof chambers or in the presence of flammable gases or fumes
- ▶ The pump must be switched OFF when changing the pump-head.

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### Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

**ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.**

ISMATEC® does not admit responsibility for the handling of chemicals.

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## Safety precautions

- ▶ The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when routing tubing in explosion-proof chambers.
- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- ▶ Do not manipulate the pump-head before the pump is switched OFF and disconnected from the mains.
- ▶ Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.

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## Safety precautions

- ▶ It is not recommended that pumped liquid contain particles. Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- ▶ Only new fuses, according to the specifications stated on Page 10 in this manual, must be used.
- ▶ The fuse-holder must not be shortcircuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies, ISMATEC® denies any responsibility.

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## Warranty

### From date of delivery

**REGLO-Z/-ZS drive:** 2 years  
**Pump-head:** 1 year  
(without consumables)

In case of any queries, please contact your nearest ISMATEC® representative.

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## Warranty Terms

We warrant the perfect functioning of the REGLO-Z/-ZS drive, provided it has been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

The pump-heads are covered by the specific terms of Micropump Ltd. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs for instruments under warranty are charged to the customer.

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## Warranty Terms

### Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use,
- ▶ unauthorized modification or misuse by the user or by a third-party,
- ▶ improper site preparation and maintenance,
- ▶ operation outside environmental and electrical specifications for the product,
- ▶ use of third-party software hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- ▶ damages caused by contamination or leaks due to torn or bursted tubing.

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## Product

### Package contents

- ▶ Gear pump drive (REGLO-Z) or control unit with gear pump drive and connecting cable (REGLO-ZS) magnet included
- ▶ Power Cord
- ▶ Operating Manual

Accessories as ordered

1. Pump-head
2. Tubing adapters

Please check the packaging and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
#### Please note

The REGLO-Z drive runs only in clock-wise direction.

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## Operating panel

1. Mains switch and Start / Stop
2. Speed selector  
99 = 5000 rpm (max. revolution)
3. Overload indicator (red LED)

## Overcurrent protector

When an overload condition occurs, an activated overcurrent protector is indicated by the red LED (3) and the pump is stopped. Whenever this situation occurs, **the pump must be switched OFF immediately**.

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again) Before the pump is re-started, it is important to check the reason for the overload (eg. too high differential pressure).

➔ Only after the cause of the overload has been detected and the failure corrected accordingly, may the pump be started again.

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## Rear panel

1. Analog interface
2. Mains socket
3. Fuse-holder with voltage selector 115/230V

## ⚠ Mains voltage

| Mains voltage           | Fuse rating     |
|-------------------------|-----------------|
| 220–240 V <sub>AC</sub> | 2 x T800mA/250V |
| 100–120 V <sub>AC</sub> | 2 x T800mA/250V |

## ⚠ Socket/Power cord

Use exclusively the originally supplied power cord. The socket must be connected to earth ground (protective conductor contact).

## ⚠ Before starting-up

Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
## Voltage setting 115V / 230V and changing the fuses

- ➔ Switch the pump OFF, pull out the  
mains plug.
1. Pull out the fuse-holder by opening it with  
a small screw-driver (size 0).
  2. Take out the voltage selector plate. Turn it  
and re-insert it into the fuse-holder so that  
the required voltage rating is facing the  
window of the fuse-holder.
  3. Insert 2 new fuses  
230 V<sub>AC</sub>: 2 x T800mA/250V  
115 V<sub>AC</sub>: 2 x T800mA/250V  
⚠ Always use 2 slow-blow fuses of the  
same type complying with the local  
mains voltage.
  4. Shut the fuse-holder. The voltage rating is  
visible in the window.

### ⚠ Before starting-up

Check that the voltage setting visible in  
the window of the fuse-holder complies  
with the local mains voltage.

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
## Please Note

Before initial start-up, the pump-head must be filled with the medium or water. Never run the pump-head under dry conditions. **Never pump liquids that contain particles. If necessary, filter the medium before use.**

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## Starting the pump

- ▶ Check the voltage rating in the window of the fuse-holder (on rear panel). If necessary, change the rating as indicated on Page 11.
- ▶ If necessary, mount the pump-head as stated on Pages 19–20.
- ▶ Before screwing the tubing adapters into the pump-head, it is recommended that **PTFE tape be wrapped around** the thread of the adapter.
- ▶ Connect the Tubing to the system. Fasten the tubing to the adapters by means of the tubing clamps.
- ▶ Connect the pump to the mains.
- ▶ Fill the pump-head with liquid before pumping starts.

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## Starting the pump

1. Set the required speed on the 2-digit speed selector. (99 = 5000 rpm)

➔ The rotation speed can be adjusted while the pump is running.

2. Press in the mains switch  
= run/stop function

➔ The drive runs only in the clock-wise direction.

⚠ **Never run the pump-head under dry conditions.**

⚠ **Never pump liquids that contain particles. If necessary, filter the medium before use.**

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
## Analog interface

### Pin 1, GND (ground)

Reference potential for all other inputs.

### Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, Pin 2 must be connected with Pin 1 (GND).

### Pin 3, start

In remote operation (Pin 2 to GND) the pump starts when connected to Pin 1 (GND).

### Pin 5, speed IN

For external speed control (0–5V<sub>ccr</sub>, 0–10V<sub>ccr</sub>, 0–20mA, 4–20mA) Input impedance and alternatives by means of DIP switches inside the pump (see Page 16).

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\*Default-Einstellung

## Settings of switch S1

➔ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected for some time.

| DIP-Switch 1 | DIP-Switch 2 | DIP-Switch 3 |
|--------------|--------------|--------------|
| ON*          | OFF*         | OFF*         |
| OFF          | ON           | OFF          |
| OFF          | OFF          | ON           |
| OFF          | OFF          | OFF          |

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| DIP-Switch 4 | DIP-Switch 5 | DIP-Switch 6 |
|--------------|--------------|--------------|
| OFF*         | OFF*         | ON*          |
| OFF          | OFF          | ON           |
| OFF          | ON           | ON           |
| ON           | ON           | OFF          |

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---

## Before switching the pump off

If a non-return valve is not mounted on the discharge side of the pump, a siphoning may occur at the moment the pump drive is stopped! (see Page 28)

If the pump will not be used for some time, it is strongly recommended that the pump-head be thoroughly rinsed with water.

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## Pumping against pressure

The max. differential pressure depends on the mounted pump-head model (see Pages 24–25).

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## Pump-head

In general, the REGLO-Z gear pump is supplied as a complete unit, consisting of the REGLO-Z drive, and as ordered:

1. pump-head
2. driving magnet
3. tubing adapters (2 pcs.)

For the REGLO-Z, a wide range of different gear pump-heads with PTFE, PPS or graphite gears are available.

An overview chart of the pump-heads available can be found on Pages 24-25.

**⚠ Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.**

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## Mounting the pump-head

### Fig. 3

Screw the 2 tubing adapters into the pump-head.  
 ➔ Before screwing the tubing adapters into the pump-head, it is recommended that **PTFE tape be wrapped around** the thread of the adapter.

### Fig. 4

Screw the pump-head to the fixing-case of the drive. The face of the pump-head with the engraved model number must be directed to the user and the arrow must point to the right.

➔ Please note that gear pumps can only pump in one direction. An arrow engraved on the front of the pump-head indicates the correct flow direction.

⚠ **Fill the pump-head with liquid before you start pumping. Avoid dry-running of the pump-head.**

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### Fig. 3

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- **»Cavity style« pump-head**
- The »Cavity style« pump-head can be used for viscous media and for applications with a certain suction lift.
- Excellent chemical resistance
- Smooth operation at a low noise level
- Low inner friction
- Incl. B mount (flange)
  
- **»Suction shoe« pump-heads**
- Enhanced pump performance at elevated differential pressures
- Fewer parts, easy to service
- Suitable for elevated temperature range
- Not suitable for applications requiring a suction lift
- Incl. B mount (flange)

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## Pump-head specifications

- **Stainless steel pump-head housing**  
316 SS/DIN 1.4571
- **Temperature range**
  - 46°C to +54°C / -51°F to +129°F (Cavity style)
  - 46°C to +177°C / -51°F to +350°F (Suction shoe) (higher temperatures possible on request)
- **Differential pressure**  
The difference between the pressure on the inlet and the discharge port of the pump-head max. 5.2 bar (75 psi) depending on pump-head.

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## Pump-head specifications

- **System pressure**  
The pressure within the pump-head  
max. 21 bar / 305 psi  
(Exception: Z-1830, Z-150: 22 bar)
- **Gear material**  
PTFE = Polytetrafluoroethylene  
PPS = Polyphenylsulfide  
Graphite = Carbon graphite  
PEEK = Polyetheretherketone
- **Ports**  
(Inside thread) 1/8" -27 NPT for  
tubing with 6.4 mm ID

The pump-heads are interchangeable. Different drives take different magnets. The REGLO-Z Analog drive is originally supplied with the correct magnet. When ordering a new pump-head, please state the intended operation temperature. If you intend to use a Cavity style pump-head for temperatures exceeding 54°C (130°F), the correct sealing material must be ordered. Please contact your nearest ISMATEC® representative. Some pump-heads are available with an **internal bypass**. This bypass is to be used for safety reasons only and **must not** be used as pressure control.

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## Suction Shoe

Für korrosive Medien  
For corrosive media  
Pour milieux corrosifs

| Pumpenkopf Nr.<br>Pump-head No.<br>Tête de pompe No.                            | Z-186       | Z-186 P      | Z-181        | Z-183        | Z-1830       | Z-1830 P     | Z-186 HC           | Z-183 HC           |
|---|-------------|--------------|--------------|--------------|--------------|--------------|--------------------|--------------------|
| Fließrate / Flow rate / Débit (ml/min) max. <sup>1)</sup>                       | 85          | 85           | 210          | 420          | 460          | 460          | 85                 | 420                |
| Fließrate / Flow rate / Débit (ml/min) min. <sup>1)</sup>                       | 0.85        | 0.85         | 2.1          | 4.2          | 4.6          | 4.6          | 0.85               | 4.2                |
| Differenzdruck<br>Pression différentielle<br>Differential pressure              |             |              |              |              |              |              |                    |                    |
|   |             |              |              |              |              |              |                    |                    |
|   |             |              |              |              |              |              |                    |                    |
| Zahnradmaterial<br>Gear material<br>Matériaux engrenage                         | Graphite    | PEEK         | Graphite     | Graphite     | PPS          | PEEK         | Graphite           | Graphite           |
| Dichtungen<br>Seals<br>Étanchéités  | PTFE        | PTFE         | PTFE         | PTFE         | PTFE         | PTFE         | PTFE               | PTFE               |
| Gehäusematerial<br>Housing material<br>Matériaux caisse                         | 316SS       | 316SS        | 316SS        | 316SS        | 316SS        | 316SS        | Hastelloy-<br>C276 | Hastelloy-<br>C276 |
| Systemdruck<br>System pressure<br>Pression de système                           |             |              |              |              |              |              |                    |                    |
|   |             |              |              |              |              |              |                    |                    |
| Betriebstemperatur<br>Operating temperature<br>Temperature operation            |             |              |              |              |              |              |                    |                    |
|   |             |              |              |              |              |              |                    |                    |
| Bypass<br>Internal bypass<br>Bypass interne                                     | --          | --           | --           | --           | --           | --           | --                 | --                 |
| Anschlüsse (Innengewinde)<br>Ports (internal thread)<br>Ports (internal thread) | 1/8" - 7NPT | 1/8" - 27NPT | 1/8" - 27NPT | 1/8" - 27NPT | 1/8" - 27NPT | 1/8" - 27NPT | 1/8" - 27NPT       | 1/8" - 27NPT       |
| Bestell-Nr.<br>Order No.<br>No. de commande                                     | MI0006      | MI0312       | MI0007       | MI0008       | MI0131       | MI0280       | MI0309             | MI0310             |

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
## Cavity Style

| Cavity Style  | Für korrosive Medien<br>For corrosive media<br>Pour milieux corrosifs |              |              |              |              |                |                |
|---|---|--------------|--------------|--------------|--------------|----------------|----------------|
|   | Z-120   | Z-140        | Z-140 P      | Z-130        | Z-150        | Z-140 HC       | Z-150 WI       |
| Pumpenkopf Nr.<br>Pump-head No.<br>Tête de pompe No.                            |   |              |              |              |              |                |                |
| Fließrate / Flow rate / Débit (ml/min) max. <sup>1)</sup>                       | 3200  | 3200         | 3200         | 3200         | 3200         | 3200           | 3200           |
| Fließrate / Flow rate / Débit (ml/min) min. <sup>1)</sup>                       | 32  | 32           | 32           | 32           | 32           | 32             | 32             |
| Differenzdruck<br>Pression différentielle<br>Differential pressure              | max. (bar)  | 3)           | 3)           | 3)           | 3)           | 3)             | 3)             |
| Zahnradmaterial<br>Gear material<br>Matériaux engrenage                         | PTFE  | PTFE         | PEEK         | PPS          | PPS          | PTFE           | NiC            |
| Dichtungen<br>Seals<br>Etanchéités  | PTFE  | PTFE         | PTFE         | PTFE         | PTFE         | PTFE           | PTFE           |
| Gehäusematerial<br>Housing material<br>Matériaux caisse                         | 316SS   | 316SS        | 316SS        | 316SS        | 316SS        | Hastelloy-C276 | Hastelloy-C276 |
| Systemdruck<br>System pressure<br>Pression de système                           | max. (bar)  | 21           | 21           | 21           | 22           | 21             | 21             |
| Betriebstemperatur<br>Operating temperature<br>Temperature operation            | °C  | -46...+54    | -46...+54    | -46...+54    | -46...+54    | -46...+54      | -46...+54      |
| Bypass<br>Internal bypass<br>Bypass interne                                     | ✓   | --           | --           | ✓            | --           | --             | --             |
| Anschlüsse (Innengewinde)<br>Ports (internal thread)<br>Ports (internal thread) | 1/8" - 27NPT  | 1/8" - 27NPT | 1/8" - 27NPT | 1/8" - 27NPT | 1/8" - 27NPT | 1/8" - 27NPT   | 1/8" - 27NPT   |
| Bestell-Nr.<br>Order No.<br>No. de commande                                     | MI0013  | MI0016       | MI0313       | MI0019       | MI0020       | MI0284         | MI0265         |

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### Operating range for cavity style pump-heads

on the REGLO-Z Analog drive: see diagram For higher pressures use the MCP-Z drive!

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- 1) Flow rate without differential pressure
- 2) For applications with differential pressures exceeding 1.0 bar (14.5 psi) we recommend the MCP-drive

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**Pump-head with bypass**

Some gear pump-head models have an internal fluid bypass.

- A fluid bypass valve helps protect against decoupling and damage from high pressure build-up.
- It allows adjustment of the max. differential pressure (from 0.7 bar up to the max. pressure valid for the individual pump-head).
- It should only be used for safety reasons and not for pressure control (bypass conditions create a temperature rise causing swelling in PTFE gears).

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## Accessories

### Service Kits:

contain the moving parts (bushes, seals, gears)

### Kit No.

41

42

82

35

36

34

88

89

## Accessories

### Tubing Adapter 01 for all gear pump-heads

1/8" thread NPT tubing I.D. 6–7 mm (1/4") (always order 2 pcs.)

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### Footswitch

Order No. ISM891

This footswitch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

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## Please Note

For maintenance, service, repairs and warranty of the pump-head we also refer to the enclosed instruction manual of the pump-head manufacturer. In case of doubt that manual shall prevail.

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## Maintenance

It is strongly recommended that the pump-head be thoroughly rinsed with water after use.

Provided the REGLO-Z/-ZS drive is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear (see Page 27).

## Repairs

For repairs, please send the defective REGLO-Z/-ZS drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., model.

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## Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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## Spare parts

If required your ISMATEC® representative will be pleased to send you on request:

- ▶ spare parts
- ▶ parts lists
- ▶ wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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## Please Note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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## Technical Specifications

### Drive

|                                     |                   |
|-------------------------------------|-------------------|
| Motor type                          | DC-Motor          |
| Speed range adjustable in 1 % steps | 50–5000 rpm       |
| Remote control                      | (see Pages 14–16) |

### Mains connection

#### Fuse rating

|             |                 |
|-------------|-----------------|
| 230 V 50 Hz | 2 x T800mA/250V |
| 115 V 60 Hz | 2 x T800mA/250V |

**Power consumption** max. 50 W

### Operating conditions

|  |               |
|--|---------------|
| Temperature  | 41 to + 104°F |
| Rel. humidity                                      | max. 80%      |
| - not condensing, at normal laboratory conditions. |               |

### Differential pressure

according to pump-head

### Dimensions/Weight

|               |                       |
|---------------|-----------------------|
|               | <b>DxWxH (inches)</b> |
| Pump drive    | 6.9x3.9x5.3           |
| Control unit  | 6.4x2.6x3.2           |
|               | <b>Weight (lb)</b>    |
| Pump drive:   | 4.4                   |
| Control unit: | 3.5                   |

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**ISMATEC®**

Learn more  
about the  
Ismatec® product  
line by visiting:  
[www.ismatec.com](http://www.ismatec.com)

**For ordering and technical support, please contact:**

**North America** [sales@ismatec.com](mailto:sales@ismatec.com) | 1-800-323-4340 | 1-847-549-7600

**Europe** [sales.europe@ismatec.com](mailto:sales.europe@ismatec.com) | +49 (0) 9377 9203-0

 **Cole-Parmer®**