

# Miniaturized Asymmetrical Flow Field-Flow Fractionation for the Fast Separation of Nanoparticle Mixtures

Dr. Florian Meier<sup>(a)</sup>, Dr. David Müller<sup>(b)</sup>, Zhengchao You<sup>(c)</sup>, Dr. Steffen Weidner<sup>(c)</sup>, Dr. Thorsten Klein<sup>(a)</sup>

<sup>(a)</sup>Postnova Analytics GmbH, <sup>(b)</sup>Centre Suisse d'Electronique et de Microtechnique (CSEM), <sup>(c)</sup>Bundesanstalt für Materialforschung und -prüfung (BAM)

## General Information

ID0049

<b>Application</b>	Nano
<b>Technology</b>	mAF4-UV
<b>Info</b>	Postnova AF2000 MT Equipped With a Miniaturized Fractionation Cartridge, PN3211 UV-Vis
<b>Keywords</b>	Miniaturized AF4, Fast Analysis Times, Nanoparticle Analysis

## Introduction

Asymmetrical Flow Field-Flow Fractionation (AF4) is a powerful tool for the separation of particles in the size range of approximately 1 - 3000 nm. However, particularly for polydisperse samples, analysis times may easily reach one hour or longer, hampering its application to high-throughput analyses. One possibility to significantly reduce analysis times in AF4 is the application of a miniaturized fractionation cartridge [1,2].

## Miniaturized Asymmetrical Flow Field-Flow Fractionation

In miniaturized AF4 (mAF4), both external dimensions as well as the separation area are significantly reduced in comparison to a conventional analytical fractionation cartridge (Figure 1).

The benefits of mAF4 include:

- Significantly reduced sample consumption and sample dilution giving improved limits of detection
- Improved compatibility with mass spectrometric detection systems due to lower flow rates
- Significantly reduced analysis time


The reduced analysis time is a major advantage of mAF4 as demonstrated in figure 2 where the elution times for a mixture of three gold nanoparticle size standards (10 nm, 30 nm and 60 nm) could be dramatically reduced in comparison to conventional AF4.

## Conclusion

Miniaturized Asymmetrical Flow Field-Flow Fractionation is a powerful complement to AF4 performed with a conventional analytical cartridge. Its low operational costs along with its fast analysis times and low sample consumption renders mAF4 a valuable tool for fast screening analyses and studies, when fully optimized peak resolution is not required.

## References

- [1] Müller D., Cattaneo S., Meier F., Welz R., deMello A.J., *Frontiers in Chemistry*, 2015, 3.
- [2] You Z., Meier F., Weidner S., *Separations*, 2017, 4(1), 8-19.



	mAF4	AF4
Dimensions [cm]	8 x 4	33.5 x 6
Separation area [cm <sup>2</sup> ]	4.3	32.8

Figure 1: Comparison of mAF4 and AF4 (with permission from [1]).




Figure 2: Comparison of fractograms obtained from AF4 (black) and mAF4 (red) for a mixture of gold nanoparticles of 10, 30 and 60 nm (with permission from [1]).