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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Caution

Rotating pumphead creates a pinch and crush hazard.

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Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and the following recommendations will be observed:

- ▶ The circuit between mains supply and pump has to be connected to earth ground.
- ▶ The pump must not be operated outside the designed operating and environmental conditions.
- ▶ The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- ▶ The pump must be switched OFF when changing the pump-head.
- ▶ Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- ▶ The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when routing tubing in explosionproof chambers.

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Please note

Never pump media containing particles.

If necessary, filter the media before starting your application.

Never run the pump-head under dry conditions.

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Safety precautions

- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- ▶ We advise you never to pump liquid containing particles. Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- ▶ Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.
- ▶ Only new fuses, according to the specifications stated on Pages 8 and 55 in this manual, must be used.
- ▶ The fuse holders must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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From date of delivery

Pump drive: 2 years
Pump-head: 1 year

The pump-heads are covered by the specific terms of Micropump Ltd. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use
- ▶ unauthorized modification or misuse by the user or by a third-party
- ▶ improper site preparation and maintenance
- ▶ operation outside of the environmental and electrical specifications for the product
- ▶ use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- ▶ MCP-Z Process drive
- 1 magnet
Order No. ISM918A
- ▶ incl. integrated power cord, 2 m long, with IEC 320 plug (male)
- ▶ 1 power cord with an IEC 320 connector (female) and a country specific mains plug
- ▶ 2 gaskets (in a plastic bag) for the analog (AD0069) and RS232-IN (AD0071) interfaces
- ▶ 2 tubing clamps
- ▶ 1 non-return valve (PP)
- ▶ Operating Manual

If ordered:

- ▶ Pump-head with 2 tubing adapters
- ▶ Software
- ▶ other accessories

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

➡ Complaints can only be accepted within 8 days from receipt of the goods.

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Please note

For using the analog or RS232 interface under IP-65 conditions please refer to Page 54.

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Socket/Power cord

The socket must be earthed. (protective conductor contact)

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Rear panel

1. Mains switch (ON/OFF)
2. Power cord
3. RS232 IN (female)
4. RS232 OUT (male)
5. Analog interface
 - ▶ input for:
 - speed control
0 – 5 V_{AC} or 0 – 10 V_{ACr} and
0 – 20 mA or 4 – 20 mA
 - Run/Stop
 - Footswitch
 - 2 digital inputs (TTL-level)
 - ▶ output for:
 - speed 0 – 10 V_{DC} or 0/12 kHz –
 - 2 universal outputs e.g. for valve unit (Pages 30 and 34)

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Fuse rating on power supply*

T4A/250V

* See also Pages 53 – 55, fuse
not present on newer models.

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Starting the pump

- ▶ Mount the pump-head as stated on Pages 45 – 46.
- ▶ Before screwing the tubing adapters into the pump-head, we recommend you wrap **PTFE tape around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.
- ▶ Connect the tubing to your system. Fasten the tubing to the adapters by means of the tubing clamps.
- ▶ Connect the pump to the mains and switch it ON with the power supply switch.

➡ Enter the ID Codes of the mounted pump-heads in the program memory currently used (see basic settings, Page 14).

⚠ **Socket/Power cord** Use exclusively the original power cord supplied with the pump. The **socket must be connected to earth ground** (protective conductor contact).

⚠ **Before you start pumping**, fill the pump-head with the medium or water. Never run the pump-head under dry conditions.

⚠ **Never pump liquids that contain particles**. If necessary, filter the medium before use.

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
4 Pro4

5

PUMP	* rpm	o PROGRAM
DISP	o Time	o Flow Rate
PAUSE	o Time	o Volume
		o TOTAL

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Before starting the pump for the first time or after connecting a new pump-head, the ID-Code must be entered in each program used with this pump-head (see basic settings, Page 14).

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Start-up information


After switching on the power supply switch, the following values are displayed:

1. LED test »8.8.8.8.«
2. Firmware version
3. Identification Code (ID-Code) of the pump-head entered in the program last used (e.g. 120)
4. Program last used
5. Displays the operating mode last used

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	○ PROGRAM
PUMP ○ rpm	● Flow rate
DISP ○ Time	● Volume
PAUSE ○ Time	● TOTAL

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
Operating panel

1. Digital LED display
2. Reduce value
* settings
3. Increase value
4. RUN/STOP
 - Starts and stops the pump
 - Interrupts and continues a dispensing cycle (Page 22)
 - ➡ The function Start/Stop can also be triggered by an external device such as a footswitch (accessories on Page 51).
5. CAL
Calibrating key for the flow rate or dispensing volume (Pages 21 and 24)
* ok
6. RESET
 - Resets a flow rate or volume calibration to the default value (Page 25)
 - Resets the accumulated volume in mode TOTAL to zero (Page 18)
 - Stops an interrupted dispensing cycle for good (Page 22)
* default
7. MAX
max. speed in the modes PUMP●rpm and ●Flow rate (keep the MAX key pressed down while the pump is running)
* escape

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Operating panel

8. MODE
Changes between the operating modes

9. LEDs for active operating mode

- ▶ **PROGRAM**
Program selection 1 – 4
- ▶ **PUMP rpm**
Pumping by speed 60 – 6000 rpm
- ▶ **PUMP Flow rate**
Pumping by flow rate in ml/min or l/min (depending on the pump-head)
- ▶ **DISP Time**
Dispensing by time 0.10s to 999h
- ▶ **DISP Volume**
Dispensing by volume in ml or litres
- ▶ **PAUSE Time**
Pause time for intermittent dispensing
- ▶ **TOTAL**
Read-out of totally delivered volume

➡ The drive runs only in the clockwise direction, there is no key for reverse rotation direction.

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Control keys for basic settings

- settings**
Allows access to the menu settings
- ok**
Confirms a data entry
- escape**
 - Cancels an entry
 - Quits the menu settings
- default**
Sets a value to its default value:
 - Cycles (dispensing cycles)
 - Head (pump-head)
 - FS (footswitch)
 - Valve
 - Digits (display)
 - Address
 - Baud Rate
- ▲▼ keys**
 - Select the basic settings
 - Select-value

➡ For selecting the basic settings refer to the following pages.

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Program specific basic settings

➔ **Individually** adjustable for each of the 4 programs.

Cycles

Number of dispensing steps by time or volume (Pages 22/23), 0..9999 (0 = Default value, infinity).

Head

Identification code (ID) of the pump head, engraved on the front of the head, see also table on Pages 48/49 (120 = default value).

➔ In order to obtain correct pumping and dispensing values, the individual ID code of the mounted pump-head must be entered; especially when a pump is used for the first time and each time the pump-head is changed.

➔ In each of the 4 programs a different pump-head can be stored.

Footswitch or hand dispenser

▶ Toggle (default value)

- Step on footswitch = pump starts running
- Press footswitch again = pump stops

▶ Direct

- the pump runs as long as the footswitch is pressed.

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Program specific basic settings

⇒ **Individually adjustable** for each of the 4 programs.

Valve (Page 30)

- ▶ Inactive (default value)
= valve is not activated
- ▶ Active
= valve is activated

Digits (display)

Number of significant digits, 3 or 4 (3 = default value)

Initializing

Pressing the ok key resets the following basic settings in the **currently used** program to the default values:

Mode:	PUMP rpm
Speed:	1000 rpm
Dispensing time:	4.50 sec.
Dispensing volume and flow rate:	depends on the pump-head
Pause time:	2.00 sec.
Number of dispensing cycles:	0

⇒ The **pump-head ID code** entered in the **currently used** program remains stored.

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Selecting the basic settings

1. Select the program (Page 18)
2. Switch the pump OFF
3. Keep the settings key pressed and switch the pump ON
4. Select the basic settings by using the ▲▼ keys

Individually adjustable for each program

- ▶ Cycles (number of disp.cycles)
- ▶ Head = ID code (pump-head)
- ▶ FS (footswitch)
- ▶ Valve
- ▶ Digits (display)
- ▶ Init (resets all basic settings to the default values)

Valid for all 4 programs

- ▶ Init4 in all 4 programs the above mentioned basic settings as well as
 - ▶ Address (RS232, pump address) and
 - ▶ Baud Rate (RS232) are reset to their default values
5. ok = confirm your selection
 6. Select the value by using the ▲▼ key, or
 - press escape for cancelling wrong input
 - certain values can be reset to the default value by using the default key
 7. ok = confirm entered value
 8. 8 escape Quit menu basic settings

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
Selecting the program

When switching the pump ON, it always selects the previously used program.

- Change mode to PROGRAM by using the MODE key.
- Press the CAL button. The currently set program starts to blink. Change to the required program with the **▲▼** keys.
- Confirm with the CAL key.
- The pump automatically returns to the last entered operating parameters of the selected program.

⚠ From now on, any changes carried out in the operating modes are automatically stored in the currently selected program.

TOTAL volume

Change with the MODE key to TOTAL. The totally delivered volume is displayed.

Four display readings are available:

Microlitre: e.g. 17.0 µl = **17.0 µ**
 Millilitre: e.g. 1700 ml = **17.00**
 Litre: e.g. 17.0 l = **17.0L**
 Litre: e.g. 1700 l = **1700.**
 (Display in litres from 9999 ml)

For zero-setting, either

- ▶ press the reset key twice, or
- ▶ switch the pump OFF.

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

- Change mode to PUMP rpm (60–6000 rpm) digitally adjustable in steps of 1 rpm
- Enter the required speed with the ▲▼ keys
Start the pump with RUN/STOP
 - ➔ The speed can also be adjusted while the pump is running.
 - ➔ During the pumping process the following values can be retrieved with the MODE-key:
 - Flow rate (ml/min)
 - Totally pumped volume (TOTAL)
- Maintain pressure on the MAX key while the pump is running.
 - ➔ Pump runs at maximum speed (ideal for fast filling or emptying the system).

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Pumping by flow rate

- Change the mode to PUMP Flow rate
 - Enter the required flow rate with ▲▼ keys (displayed in ml/min)

Start the pump with RUN/STOP


 - ➔ The flow rate can also be adjusted while the pump is running.
 - ➔ For an accurate flow rate it is recommended that the pump should be calibrated (see Page 21).
 - ➔ During the pumping process the following values can be retrieved with the MODE key:
 - Rotation speed (PUMP rpm)
 - Totally pumped volume (TOTAL)
- Maintain pressure on the MAX key while the pump is running
 - ➔ Pump runs at maximum speed (for fast filling or emptying the system)

⚠ Entering the ID-Code of the currently mounted pump-head (Pages 48/49) allows you set the flow rate in ml/min. This, however, is an approximate value as not yet calibrated.

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Calibrating the flow rate

- Change the mode to PUMP Flow rate. Enter the required flow rate with the ▲▼ keys.
 - Enter 60 seconds by using the ▲▼ keys
 - Start the pump with RUN/STOP

➡ The pump stops automatically after 60 sec.

➡ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:
- Change the mode to DISP Time.
- Return to the mode PUMP Flow rate.
- Press the CAL button (displayed value blinks) Enter the weighed or measured value with the ▲▼ keys.
- Confirm with the CAL key (the flow rate setting returns automatically to the initially entered set point).
- Start with the RUN/STOP key.
 - ➡ Depending on the application conditions repeated calibration may be necessary.

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Please note

With RUN/STOP a dispensing cycle can be interrupted and resumed. With RESET an interrupted dispensing cycle can be definitely stopped.

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Dispensing by time

The dispensing time can be entered from 0.10 s – 999 h.

1. Change mode to DISP Time.
2. Enter the required dispensing time with the ▲▼ keys.

➡ The display accelerates when pressure on the ▲ or ▼ key is maintained.

Seconds: 0.10" – 99.9" (in 0.01 s steps)
0.1" – 899.9" (in 0.1 s steps)
Minutes: 15' – 899' (in 1 min steps)
Hours: 15 h – 999 h (in 1 h steps)

3. Start pump with RUN/STOP.

➡ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point only takes effect at the beginning of the subsequent dispensing cycle.

➡ During the dispensing process the following values can be retrieved by pressing the MODE-key:

- speed (PUMP rpm)
- flow rate (PUMP Flow rate)
- totally dispensed volume (TOTAL)

➡ With the ▲▼ keys the speed and flow rate can be changed even during the dispensing process.

➡ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible.

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PUMP	o rpm	o PROGRAM
DISP	o Time	o Flow Rate
PAUSE	o Time	* Volume
	o TOTAL	

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Dispensing by volume

1. Change the mode to DISP Volume.

2. Use the ▲▼ keys for entering the required dispensing volume. Four display readings are available:

Microlitre: e.g. 17.0 µl = **17.0µ**
Millilitre: e.g. 1700 ml = **17.00**
Litre: e.g. 17.0 l = **17.0L**
Litre: e.g. 1700 l = **1700**.
(Display in litres from 9999 ml)

➔ For dispensing accurately, it is recommended that the pump should be calibrated (Page 24).

3. Start the pump with RUN/ STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

➔ The dispensing speed can be adjusted in the modes PUMP rpm and PUMP flow rate.

➔ With the ▲▼ keys the dispensing volume can be changed even during the dispensing process. The newly entered set point appears shortly in the display.

➔ During the dispensing process the following values can be retrieved with the MODE key:

- Rotation speed (PUMP rpm)
- Flow rate (PUMP Flow rate)
- totally dispensed volume (TOTAL)

➔ With the ▲▼ keys the rotation speed or flow rate can be changed even during the dispensing process.

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Calibrating the volume

1. Change the MODE to DISP Volume.

➔ Use the ▲▼ keys for entering the required dispensing volume.

2. Start the pump with RUN/STOP.

➔ Pump stops automatically.

➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes.

➔ Measure the dispensed liquid by volume or weight.

3. Stay in the same mode (DISP Volume) and press the CAL-key (display blinks).

4. Enter the effectively dispensed volume with the ▲▼ keys.

5. Save the entered value with the CAL key.

➔ **Pump is calibrated and returns automatically to the initially entered set point.**

6. Start the pump with RUN/STOP.

➔ Depending on the application conditions repeated calibration may be necessary.

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1 **MODE**

PUMP	rpm	PROGRAM
DISP	Time	*Flow Rate
PAUSE	Time	o TOTAL

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Default calibration of flow rate

1. Change the mode to PUMP Flow rate.
 2. Press the CAL button (displayed value blinks).
 3. Press the RESET button (displayed value blinks).
 4. Confirm with the CAL button.
- ➡ At the same time, the volume is also reset to the default value.

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Default calibration of volume

1. Change the mode to DISP Volume.
 2. Press the CAL button (displayed value blinks).
 3. Press the RESET button (displayed value blinks).
 4. Confirm with the CAL button.
- ➡ At the same time, the flow rate is also reset to the default value.

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Dispensing by volume within a pre-set time

1. Change the mode to DISP Time. Enter the required time with the ▲▼ keys (0.10 s – 999 h).

2. Change the mode to DISP Volume. Enter the required volume in ml with the ▲▼ keys.

3. Push the MODE button.

➔ Both LEDs light up simultaneously.

4. Start the pump with RUN/STOP.

➔ The pump dispenses the required volume within the pre-set time.

a Changes / Corrections
Press the Reset key

➔ For changing the dispensing time or dispensing volume, the pump returns directly to the mode DISP Time

➔ Repeat steps 1 to 3

b Error message

LEDs are blinking alternately

9999 = volume too large time too short
1111 = volume too small time too long

Despite volume and time adjustments, the LEDs keep blinking:

➔ Change the pump-head model in order to obtain another flow rate range.

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Please note

With RUN/STOP a dispensing or pause cycle can be interrupted and resumed.

With RESET an interrupted dispensing or pause cycle can be definitely stopped.

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Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time.

- Go to the mode PAUSE Time. Enter the pause time with the ▲▼ keys (between 0.10s – 999h).
- Press the MODE key once
 ➔ Both LEDs for PAUSE and DISP Time light up simultaneously. Enter the dispensing time with the ▲▼ keys (0.10 s – 999 h).
- Start the pump with the RUN/STOP key.
 ➔ Pump dispenses atm intervals until the RUN/STOP button is pushed again.
 ➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.

➔ During the dispensing process the following values can be retrieved with the MODE-key:

- rotation speed (PUMP rpm)*
- flow rate (PUMP Flow rate)*
- totally dispensed volume (TOTAL)

* The rotation speed and flow rate can be changed with the ▲▼ keys.

For entering the number of dispensing cycles see Page 29.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time.

1. Change mode to PAUSE Time. Enter the pause time (between 0.10 s – 999 h) with the ▲▼ keys.

2. Press the MODE-key twice.

➔ Both LEDs for PAUSE and DISP Volume light up simultaneously. Enter the volume with the ▲▼ keys.

3. Start the pump with the RUN/STOP key.

➔ Pump dispenses at intervals until the RUN/STOP button is pushed again.

➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing volume or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.

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- speed (PUMP rpm)*
- flow rate (PUMP Flow rate)*
- totally dispensed volume (TOTAL)
- * The rotation speed and flow rate can be changed with the ▲▼ keys.

If necessary, calibrate the volume as indicated on Page 24.

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PUMP	◻ rpm	◻ PROGRAM
DISP	◻ Time	◻ Flow Rate
PAUSE	* Time	◻ Volume
		◻ TOTAL

PUMP	◻ rpm	◻ PROGRAM
DISP	◻ Time	* Volume
PAUSE	* Time	◻ TOTAL

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

See basic settings on Page 14.

➡ During the pauses the display shows the remaining number of dispensing cycles.

➡ During long pauses, the remaining interval time can be checked with the CAL key and adjusted by using the ▲▼ keys.

➡ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key. When resuming the process, the pump continues with the remaining number of dispensing cycles.

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
Dispensing with a valve

Using a valve prevents a **syphoning effect** (liquid flows back to the reservoir) at the moment the pump is stopped. The valve opens and closes simultaneously when the pump is started and stopped. As such it is essential for accurate and repeatable dispensing.

- ▶ Connect the valve to the analog interface on the rear panel of the MCP-Z Process.
- ▶ Activate the valve
 - ➔ see basic settings on Page 15.

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Pumping against pressure

The max. differential pressure depends on the mounted pump-head model (see Pages 48/49).

For our pressure control unit, see Page 50.

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Pump-head with bypass

Some gear pump-head models have an internal fluid bypass.

- ▶ A fluid bypass valve helps protect against decoupling and damage from high pressure build-up.
- ▶ It allows adjustment of the max. differential pressure (from 0.7 bar up to the max. pressure valid for the individual pump-head).
- ▶ It should only be used for safety reasons and not for pressure control (bypass conditions create a temperature rise causing swelling in PTFE gears).

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Before switching the pump off

⚠ If no non-return valve is mounted on the discharge side of the pump, a siphoning may occur at the moment the pump drive is stopped! See Page 30.

If the pump will not be used for some time, it is strongly recommended that the pump head be rinsed with tap water.

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Overcurrent protector

The drive MCP-Z *Process* features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL S = Speed set-point cannot be reached

OL I = Current Overload

Whenever this situation occurs, the pump must be switched OFF immediately. Before the pump is restarted, it is most important to check the reason for the overload (eg. too high differential pressure, etc.).

➡ Only when the cause of the overload has been detected and the failure corrected accordingly, should the pump be started again.

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
Pin 2, remote
 Pin 3, start
 Pin 6, input 1
 Pin 8, input 2
 Pin 13, speed internal

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Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, Pin 2 must be connected with Pin 1 (GND).

Pin 3, start

- ▶ In normal operation (Pin 2 open) Pin 3 operates the footswitch (settings see Page 14)
- ▶ In remote operation (Pin 2 to GND) the pump starts when connected to Pin 1 (GND)

Pin 5, speed IN

External speed control (0 – 5 V_{DC}, 0 – 10 V_{DC}, 0 – 20 mA, 4 – 20 mA) Input impedance and input range can be selected via a DIP-switch inside the pump (see Page 35)

Pin 6, input 1 Pin 8, input 2

Digital inputs (TTL-level). Active when connected to Pin 1. Are only retrievable via RS232 interface; e.g. with a software like the ProgEdit.

Pin 7, +36 V_{DC}

About +36 V_{DC} are available (max. current 1A).

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Analog interface

Pin 9, speed OUT

The default setting is 0 – 10 V_{DC}, proportionally to the motor speed 0 – 6000 rpm.

Alternatively a frequency range from 0/12 kHz is available.

Alternatives with sliding switch S2 inside the pump (see Page 35).

Pin 10, +5 V_{DC}

About +5 V_{DC} are available. (max. current 0.5 A)

Pin 13, speed internal

Depending on how the pump is operated, Pin 13 has different functions:

Analog interface not activated

(Normal operation, i.e. Pin 2 is open) Pin 13 serves as **auto-start function**. If Pin 13 is connected to Pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be ON).

Analog interface activated
 (Pin 2 on GND)

- ▶ Pin 13 open: The rotation speed is adjusted via Pin 5 (speed IN).
- ▶ Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

Pin 14, output 1 (PWM)

Pin 15, output 2 (PWM)

Digital outputs in connection with Pin 7 (+36 V_{DC}), implemented as **pulse-width-modulation** (open collector). Can only be activated via RS232 interface, e.g. with a software like the ProgEdit.

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Please note

The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected for some time.

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This switch affects Pin 9, motor speed

Position A: 0 – 10 VDC (Standard)
Position B: 0/12 kHz

➔ In order to maintain the optimal casing seal, the pump should not be opened without a particular reason. This is, however, necessary in order to gain access to the switches. For this purpose we recommend that you proceed exactly as indicated on Pages 53 – 54.

➔ **Make sure that the pump is disconnected from the mains supply.**

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Please note

When switching the MCP-Z Process drive ON, the serial interface will respond after a 3 second delay.

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Serial interface

RS232 IN (female)

A 9-Pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1 – 8).

It is factory-set at 1 and can be changed with @ (also refer to the basic settings described on Page 16).

The address will be stored permanently (even after the pump has been switched OFF). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed** by a carriage return (ASCII 13). The pump confirms most of the commands with an asterix*. Yes/No inquiries are answered by + (yes) or – (no).

Multidigit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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Cascading several pumps

Providing that appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-Pin RS232 connector.

If the PC is equipped with a 25-Pin connector, you need an additional adapter XC0009.

➡ **Each pump must be allocated its own individual address (see also the basic settings on Page 16).**

➡ Alternatively, the addressing can also be accomplished via the serial interface.

For that, each pump must first be connected individually to a PC and the address must be allocated with @. Only then may the pumps be daisy-chained.

Each pump requires a definite address in order to avoid interference of their responses.

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Programming Software

Software ProgEdit

- ▶ Create an individual application program.
- ▶ Download the program via RS232 interface to your MCP-Z Process pump drive.
- ▶ Now run your application PC-independently wherever needed.

System requirements:

3 1/2" floppy drive, Win95 or higher, Pentium or higher.

Order No. SOF104

LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments. <http://www.ni.com/>

LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our website: <http://www.ismatec.com>

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9 pins, female / 9 pins, male (2 m long)

Order No. AG0013

(Included with the delivery of the software)

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Pump-head

In general, the MCP-Z Process gear pump is supplied as a **complete unit**, consisting of the MCP-Z Process drive, and as ordered:

1. pump-head
2. driving magnet
3. tubing adapters (2 pcs.)

For the MCP-Z Process, a wide range of different gear pump-heads with PTFE, PPS or graphite gears are available.

You find an overview chart of the pump-heads available on Pages 48/49.

⚠ Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.

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Mounting the pump-head

Fig. 1

The following tools are required for mounting a pump-head:

- ▶ Allen key 3/32"
- ▶ Phillips screwdriver (Size 0)

Remove all screws from the fixing case (a) of the drive.

(Together with the drive, two long screws are additionally supplied for fastening the large pump-heads of the Series 200 at the top and the bottom of the pump-head face).

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Fig. 1

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Mounting the pump-head

Fig. 2

Fasten the driving magnet on the motor shaft (through the hole underneath).

➔ Ensure that there is a gap of about 3 mm between the magnet and the brim of the fixing case.

⚠ **During operation the back of the magnet must not touch the housing of the pump drive.**

Fig. 3

Screw the 2 tubing adapters into the pump-head.

➔ Before screwing the tubing adapters into the pump-head, it is recommended that **PTFE tape be wrapped around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.

Fig. 4

Screw the pump-head to the fixing case of the drive. The face of the pump-head with the engraved model number must be oriented toward the user. The engraved arrow indicates the flow direction.

➔ Please note that gear pumps can only pump in one direction. An arrow engraved on the front of the pump-head indicates the correct flow direction.

⚠ **Fill the pump-head with liquid before pumping begins. Avoid dry-running of the pump-head.**

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Fig. 2

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»Cavity style« pump-head

- ▶ The »Cavity style« pump-head can be used for viscous media and for applications with a certain suction lift
- ▶ Excellent chemical resistance
- ▶ Smooth operation at a low noise
- ▶ Low inner friction
- ▶ Incl. B mount (flange)

»Suction shoe« pump-heads

- ▶ Enhanced pump performance at elevated differential pressures
- ▶ Fewer parts, easy to service
- ▶ Suitable for elevated temperature range
- ▶ Not suitable for applications requiring a suction lift
- ▶ Incl. B mount (flange)

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Pump-head specifications

- ▶ Stainless steel pump-head housing SS 316/DIN 1.4571

- ▶ **Temperature range**
- 46°C to +54°C / -51°F to +129°F (Cavity style)
- 46°C to +177°C / -51°F to +350°F (Suction shoe) (higher temperatures possible on request)

- ▶ **Differential pressure**
The difference between the pressure on the inlet and the discharge port of the pump-head max. 5.2 bar (75 psi) depending on pump-head (Pages 48/49)

- ▶ **System pressure**
The pressure within the pump-head max. 21 bar / 305 psi (Exception: Z-1830, Z-150: 22 bar)

- ▶ **Gear material**
PTFE = Polytetrafluoroethylene
PPS = Polyphenylenesulphide
Graphite = Carbon graphite
PEEK = Polyetheretherketone

- ▶ **Ports**
(Inside thread) 1/8" -27 NPT for tubing with 6.4 mm ID

The pump-heads are interchangeable. Different drives take different magnets. The MCP-Z Process drive is originally supplied with the correct magnet.

When ordering a new pump-head, please state the intended operation temperature. If you intend to use a Cavity style pump-head for temperatures exceeding 54°C (130°F), the correct sealing material must be ordered. Please contact your nearest ISMATEC representative.

Some pump-heads are available with an **internal bypass**. This bypass is to be used for safety reasons only and must not be used as pressure control (see Page 31).

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Pumpenkopf ID Code Pump-head ID Code Tête de pompe ID Code	181	183	183HC	186	186P	186HC	1830	1830 P	200	200-P	200 PKC	201	201 PKC
Bestell-Nr. Order No. No. de commande	MI0007	MI0008	MI0310	MI0006	MI0312	MI0309	MI0131	MI0280	MI0022	MI0306	MI0341	MI0023	MI0378
Fließrate Flow rate Débit ml/min max.)	252	504	504	102	102	102	550	550	3380	3380	3380	7020	7020
Fließrate Flow rate Débit ml/min max.)	2.5	5.0	5.0	1.0	1.0	1.0	5.5	5.5	34	34	34	70	70
Differenzdruck max. Pression différentielle Differential pressure (bar)	2.8	2.8	2.8	1.4	2.3	1.4	5.2	5.2	5.2	5.2	5.2	3.5	3.5
Differenzdruck max. Pression différentielle Differential pressure (psi)	40.6	40.6	40.6	20.3	33.3	20.3	75.4	75.4	75.4	75.4	75.4	50.8	50.8
Zahnradmaterial Gear material Matériaux engrenage	Graphite	Graphite	Graphite	Graphite	PEEK	Graphite	PPS	PEEK	PPS	PEEK	PEEK	PPS	PEEK
Dichtungen Seals Etanchéités	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	Viton	Viton	Kalrez® **)	Viton	Kalrez® **)
Gehäusematerial Housing material Matériaux caisse	SS316	SS316	*) Hastelloy C276	SS316	SS316	*) Hastelloy C276	SS316	SS316	SS316	SS316	SS316	SS316	SS316
Systemdruck max. System pressure Pression de système (bar)	21	21	21	21	21	21	21	21	21	21	21	21	21
Systemdruck max. System pressure Pression de système (psi)	305	305	305	305	305	305	305	305	305	305	305	305	305
Betriebstemperatur Operating temperature Temperature operation (°C)	-46...+177	-46...+177	-46...+177	-46...+177	-46...+177	-46...+177	-46...+177	-46...+177	-29...+177	-29...+177	-29...+177	-29...+177	-29...+177
Betriebstemperatur Operating temperature Temperature operation (°F)	-50...+350	-50...+350	-50...+350	-50...+350	-50...+350	-50...+350	-50...+350	-50...+350	-20.2...+350	-20.2...+350	-20.2...+350	-20.2...+350	-20.2...+350
Bypass Internal bypass Bypass interne	--	--	--	--	--	--	--	--	▼	--	▼	▼	▼
Anschlüsse (Innengewinde) Ports (internal thread) Ports (internal thread)	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT

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Pumpenkopf Nr. Pump-head No. Tête de pompe No.		120	130	Z-140	140 P	140 HC	142	142 HC	150	150 WI
Bestell-Nr. Order No. No. de commande		MI0013	MI0019	MI0016	MI0313	MI0284	MI0018	MI0311	MI0020	MI0265
Fließrate Flow rate Débit	ml/min max. ¹⁾	3840	3840	3840	3840	3840	5640	5640	3840	3840
Fließrate Flow rate Débit	ml/min max. ¹⁾	38	38	38	38	38	56	56	38	38
Differenzdruck max. Pression différentielle Differential pressure	(bar) (psi)	3.5 50.8	5.2 75.4	3.5 50.8	5.6 81.2	3.5 50.8	3.5 50.8	3.5 50.8	5.2 75.4	5.2 75.4
Zahnradmaterial Gear material Matériaux engrenage		PTFE	PPS	PTFE	PEEK	PTFE	PTFE	PTFE	PPS	Nickel-Carbid
Dichtungen Seals Étanchéités		PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE
Gehäusematerial Housing material Matériaux caisse		SS316	SS316	SS316	SS316	**) Hastelloy C276	SS316	**) Hastelloy C276	SS316	***) Oberfläche Surface hardened
Systemdruck max. System pressure Pression de système	(bar) (psi)	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6	21 304.6
Betriebstemperatur Operating temperature Temperature operation	°C °F	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129	-46...+54 -50...+129
Bypass Internal bypass Bypass interne		✓	✓	--	--	--	--	--	--	--
Anschlüsse (Innengewinde) Ports (internal thread) Ports (internal thread)		1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT

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Footswitch

Order No. IS10039

This footswitch serves as a start/stop device. Use of the footswitch allows both hands to be free to handle bottles and tubing when using this pump as a dispenser.

Connecting cable

For connecting an MCP-Z Process to a PC via the 9-Pin RS232 connector.

Order No. AG0013

Service Kit

The service kit provides all spare parts in case gears and seals require replacement.

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Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

After use, it is recommended that the pump head be flushed with warm tap water.

Provided the MCP-Z *Process* is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear (see Page 51).

Repairs

For repairs, please send the MCP-Z *Process* drive to your ISMATEC® representative. Please give information on defect, date of purchase, serial No., and model.

Spare parts

For repairs carried out by the owner (out of the warranty period) please ask your ISMATEC® distributor for:

- ▶ spare parts
- ▶ parts lists
- ▶ wiring diagrams

Please give information on defect, date of purchase, serial No., and model.

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Removing the casing hood

 The pump should only be opened for replacing the fuses or changing the DIP-switches. In order to maintain the best possible seal of the enclosure, the following steps are recommended:

Removing the casing hood

1. Remove the 4 screws with a Phillips screw driver (size 2) on each side of the casing hood.
2. Spread **slightly** the casing hood at the bottom and remove it by carefully pulling it upwards.

Mounting the casing hood

1. Before replacing the casing hood we recommend lubricating the sealing on the contact surface (e.g. with silicone oil spray - order No. SC0179).

In case that the two sealing ends have come apart, they must be stuck together again with Silicone sealing compound.

2. Spread the casing hood slightly at the bottom and push it carefully over the pump. Be careful that no cables are trapped.
3. Press the casing hood downwards and fasten the screws evenly from top to bottom.

When tightening the screws, be careful that the profile sealing is not trapped. If the sealing should be damaged, it must be replaced, otherwise the seal may no longer be tight.

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
Removing the casing hood

Sealing the casing hood

When the casing hood fits tight, each single screw must be removed again and stopped with silicone sealing compound. Then, the screws must be refastened again evenly from top to bottom.

➔ Remove the screw - fill the hole with silicone - fasten the screw again - remove the last bit of silicone with a rag. If necessary clean up with acetone.

Plug gaskets

If the interface ports are used under **IP-65 conditions**, the originally supplied gaskets must be inserted into the male plugs of the corresponding connecting cables.

➔ The interface »RS232 OUT« (male) is already equipped with a gasket.

1. Gasket for plug (15 pin, male) for analog interface (female)
Spare part No. AD0069
2. Gasket for plug (9 pin, male) for RS232-IN interface (female)
Spare part No. AD0071

⚠ **Choosing the correct interface cable (sealed plug, e.g. moulded plug) according to the IP-rating used, is entirely the user's responsibility.**

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

Changing the fuses

The fuses are fixed to the control board (above) and the power supply circuit board (below) as illustrated opposite.

Fuse rating on control board

- a. F1.6A/250V
- b. T4A/250V

Fuse rating power supply*

- c. T4A/250V

*** Fuse not present on newer models.**

⚠ Switch the pump OFF, pull out the mains plug.

⚠ Remove the casing hood. For removing and remounting the hood, please refer to Pages 53 – 54.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical specifications

Drive

Motor type DC-Motor
Speed range 60 – 6000 rpm
digitally adjustable in steps of 1 rpm

Differential pressure

max. 6 bar (87 psi) depending
on pump-head (Page 31)

Remote control

via RS232 and analog interface
(Pages 33 and 37)

Mains connection

100 - 240 V_{AC} / 50/60 Hz

Fuse rating on control board

F1.6A/250V
T4A/250V

Fuse rating power supply*

T4A/250V

* Fuse not present on newer models.

Power consumption

Max. 150 W

Protection rating

IP 65

Operating conditions

Temperature 41°F to 104°F
Rel. humidity max. 80% at
normal environmental conditions

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


Centering flange

The centering flange fits onto MCP-Z drive. It is only required for the Micropump 200 and 201 pump-head series where it facilitates their precise mounting.

Accessories:

- a. 4 flat-headed screws
M 2.5 x 8 BN0187
- b. 1 centering flange IS10379

Mounting instructions:

1. Please ensure that the circular drive magnet (c) is correctly mounted.
2. When mounting the centering flange, make sure that:
 - its elevated inner ring (1) exactly fits into the drive flange, and
 - the two straight holes (2) are placed in a vertical position at the top and the bottom.

➔ When screwing the centering flange onto the drive, make certain that none of the four flat-headed screws protrude.
3. Mounting the gear pump-head with the two screws can now proceed in the usual way.

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ISMATEC®

Learn more
about the
Ismatec® product
line by visiting:
www.ismatec.com

For ordering and technical support, please contact:

North America sales@ismatec.com | 1-800-323-4340 | 1-847-549-7600

Europe sales.europe@ismatec.com | +49 (0) 9377 9203-0

 **Cole-Parmer®**