

Stop CE Method Transfer Anxiety... We Can Make **Your** Buffers for You!

Rationale:

CE Methods may not “transfer” well because making a CE Buffer is not a trivial matter. The water you start with, the method of reagent addition, pH and BGE control are only some of the factors that can vary from lab to lab. Following and validating an SOP for buffers can be tedious and difficult; creating resistance to the method as well as reducing its ruggedness.

Transfer your CE Method successfully and easily by making it easier for other labs to perform your method with your buffer. Simply provide these labs with a catalog number to order your “Custom Made Buffer”. The labs will receive your buffer with “exact intended properties for CE” the first time, every time...worldwide.

How Do “Custom Buffers” work?

You supply MicroSolv with your buffer recipe, which we keep confidential. Your buffers are manufactured for you in small production batches which include Certificates of Analysis, Commercial Packaging and MSDS sheets when requested. We supply you with a discreet catalog number which you can refer to anyone for their use. Through our international sales network, these buffers are available worldwide and are exactly the same buffers you used to develop your original CE or CEC method. This eliminates one reproducibility variable from your method, making it more robust and acceptable to labs that are not expert in preparing CE buffers.

Eliminate “Buffer Stress”.

Ease method transfer challenges and improve chances of success for CE methods that you develop and export to another laboratory. Most analysts are not experienced with CE buffers and often do not make their buffers in the same way that you do. This results in reproducibility problems and can result in resistance to using your method.



- **Low Cost, Low Minimums**
- **Your Buffer Made in a GMP Facility**
- **Standard Catalog No. and SOP – Worldwide**
- **Easy, Sensible and Consistent**

Q. How do I find out how to get my buffers made?

A. Three easy ways.

1. **Internet.** Visit our website and fill out our “Custom Buffers” form.
2. **Phone.** Give us a call, we will gladly walk you through the process.
3. **Fax.** Send us an idea of what you are looking for and one of our representatives will be happy to fax back a proposal.

Your choice above will only take a few minutes to complete and there is no obligation to buy. Give us a try!